

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/015764

International filing date: 06 May 2005 (06.05.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US
Number: 60/569,370
Filing date: 06 May 2004 (06.05.2004)

Date of receipt at the International Bureau: 20 June 2005 (20.06.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1331928

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

June 09, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

APPLICATION NUMBER: 60/569,370

FILING DATE: *May 06, 2004*

RELATED PCT APPLICATION NUMBER: *PCT/US05/15764*



Certified by

Under Secretary of Commerce
for Intellectual Property
and Director of the United States
Patent and Trademark Office



16076 U.S. PTO

Express Mail Label No. EL988555347US

Please type a plus sign (+) inside this box → ☐

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Approved for use through 7/31/2003.

6-03

PROVISIONAL APPLICATION FOR PATENT COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

22151 U.S. PTO
60/569370

INVENTOR(S)					
Given Name (first and middle [if any])		Family Name or Surname		Residence (City and either State or Foreign Country)	
Trent Russell Neal Walter		Nortnen Woodbury		Tempe, Arizona Tempe, Arizona	
Additional inventors are being named on the ___ separately numbered sheets attached hereto					
TITLE OF THE INVENTION (500 characters max)					
LIGHT DIRECTED SOLID PHASE SYNTHESIS ON PATTERNED POLYMERS					
Direct all correspondence to: CORRESPONDENCE ADDRESS					
<input checked="" type="checkbox"/> Customer Number <input type="text" value="26707"/>					
OR Type Customer Number here					
<input type="checkbox"/> Firm or Individual Name					
Address					
Address					
City		State		ZIP	
Country		Telephone		Fax	
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification		Number of Pages <input type="text" value="54"/>		<input type="checkbox"/> CD(s), Number <input type="text"/>	
<input type="checkbox"/> Drawing(s)		Number of Sheets <input type="text"/>		<input type="checkbox"/> Other (specify) <input type="text" value="Cover sheet; Postcard"/>	
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76					
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT					
<input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.				FILING FEE AMOUNT (\$)	
<input type="checkbox"/> A check or money order is enclosed to cover the filing fees				<input type="text" value="\$80"/>	
<input checked="" type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: <input type="text" value="17-0055"/>					
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.					
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.					
<input checked="" type="checkbox"/> No.					
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____					

Respectfully submitted,

SIGNATURE

TYPED or PRINTED NAME

Robert D. Atkins

TELEPHONE

(602) 229-5311

Date

05/ 6 / 04

REGISTRATION NO.
(if appropriate)
Docket Number:

34,288

112624.00138


USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is used by the public to file (and by the PTO to process) a provisional application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the complete provisional application to the PTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Patent Application, Commissioner for Patents, Alexandria, VA 22313-1450.

EXPRESS MAIL CERTIFICATE: EL988555347US

I hereby certify that this correspondence listed below is being deposited with the United States Postal Service on the date set forth below as Express Mail in an envelope addressed to: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

Date of Signature
and Deposit: May 6, 2004

By: 
(Signature of person depositing mail)
MARITZA O'NEILL

CERTIFICATE OF MAILING PURSUANT TO 37 C.F.R. 1.10

Applicant: *Northen et al.*

Filed: May 6, 2004

Title: *LIGHT DIRECTED SOLID PHASE
SYNTHESIS ON PATTERNED POLYMERS*

Docket No.: 112624.00138

Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

Type of Filing:

- 1) Provisional Application For Patent Cover Sheet
- 2) Specification (54 pages, plus cover sheet)
- 3) Return postcard

This Page Is Inserted by IFW Operations
and is not a part of the Official Record

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

IMAGES ARE BEST AVAILABLE COPY.

**As rescanning documents *will not* correct images,
please do not report the images to the
Image Problem Mailbox.**

PATENT

PROVISIONAL APPLICATION

of

TRENT RUSSELL NORTHEN

NEAL WALTER WOODBURY

For

UNITED STATES LETTERS PATENT

on

LIGHT DIRECTED SOLID PHASE SYNTHESIS ON PATTERNED POLYMERS

Attorneys:

QUARLES & BRADY STREICH LANG L.L.P.

ONE RENAISSANCE SQUARE

TWO NORTH CENTRAL AVENUE

PHOENIX, AZ 85004-2391

Express Mail Label No.: EL988555347US
Attorney Docket No.: 112624.00138

(Include additional names and addresses on a separate sheet.)

II. DESCRIPTIVE TITLE OF INVENTION

Light directed solid phase synthesis on patterned polymers

III. GRANT/CONTACT (If any):

Sponsor(s) none Award Number _____

Principal Investigator: n/a ORSPA Acct Number _____

IV. LAB/DEPARTMENT WHERE DEVELOPED:

Woodbury lab, Department of Chemistry and Biochemistry, ASU

V. DESCRIPTION OF INVENTION:

A. This invention is a(n): x process _____ chemical compound

_____ electronic circuit _____ mixture of chemical compounds

_____ apparatus _____ therapeutic method _____ other

(describe)

B. State, as fully as possible, what the invention is, including: materials or components used; operative and preferred ranges of process parameters and concentrations of chemical compounds; and foreseeable uses of the invention.

This invention combines three existing technologies: 1. photopolymers (photoresist), 2. photolabile protective groups, and 3. solid phase synthesis. The combination of these three existing technologies allows for the construction of three dimensional surfaces and devices that have tailored chemical functionality in spatially defined areas. The substrate can be any polymer that is in the appropriate three dimensional form having groups that can be derivatized in a way that can then be protected with photolabile protective groups. Photolabile protective groups can

three dimensional form having groups that can be derivatized in a way that can then be protected with photolabile protective groups. Photolabile protective groups can include any group that can be removed with light or activated by light in a way to expose or react with a material introduced in solution. Ways of patterning the polymer may include photopolymerization, thermal polymerization, or contact stamping. Ways of removing the photoprotective group include using a scanning laser system, micromirror array, or photolithographic method. Compounds that can be attached to this surface can be almost anything that will react with the given functionality exposed upon removal of the photolabile protective groups. It is possible to use both single and multiphoton excitation of the polymer and protective group to generate the spatial features. Sequential steps of removing the photolabile protective group and coupling new materials with the protective group blocking the appropriate reactive groups, can be used to generate complex patterns of functionalized polymer surfaces. In the short term these would be useful for enhanced DNA and Peptide microarrays, longer term these could be used for things as diverse as drug delivery systems, sensors, and artificial organs.

This includes a vast number of materials and methods. Work to date has been with acrylates and methacrylates including ones with reactive side chains (epoxy) that can be functionalized with diamines to yield aminated surfaces. These polymer surfaces have been three dimensionally patterned from photoreactive monomer/polymer solutions using a scanning laser system on top of a methacrylate functionalized glass surface. They have been protected with the nitroveratryloxycarbonyl (NVOC) photolabile protecting group (there are several protecting groups for amines and hydroxyls based on the general class of nitrobenzene compounds, though there are other classes of protective groups that can be used). This protective group has been removed using a scanning laser system. Detection of the deprotective areas has been done using fluorescence from a reactive dye that selectively reacts with the exposed amine groups.

- C. Records Supporting Invention: Identify records which establish dates of conception and reduction to practice, including identity of person who prepared record and its present location. Attach copies if possible. Note additional supporting evidence. If the invention or a significant aspect of the invention is not supported by written records, briefly describe how the date of invention can be established and identify earliest written record.


A glass cover slide was cleaned for 15 min at RT with 60/40 sulfuric acid/hydrogen peroxide, placed in 10% sodium hydroxide at 70 C for 3 min, placed in 1% HCl at RT for 1 min, between each step it was soaked in nanopure water for 3 minutes. A solution of 1% 3-(trimethoxysilyl)propyl methacrylate in 95% ethanol 5 % water was made and mixed for 10 minutes, the slide was then added and left to react at RT for 15 minutes with gentle agitation. This slide was soaked in isopropyl alcohol for 3 min then nanopure water for 1 min then placed in a 100 C oven for 5 minutes after which the oven was

turned off and nitrogen was blown through for 1 hr. The slide was stored under nitrogen until it was used.

A blend of methacrylate monomers and photoinitiator (900 uL trimethylolpropane trimethacrylate, 100uL glycidyl methacrylate, 10 mg azobisisobutyronitrile) was prepared and nitrogen was bubbled through for 15 minutes before it was injected into a Focht cell (Bioptechs inc. Butler, PA) that had been flushed with argon. This cell was then mounted onto a Prior scientific microscope stage on a Nikon microscope and illuminated through a 40x 0.75 NA objective with 370nm light from a Ti:Saphire laser modulated by a Conoptics modulator (shutter). The system was controlled via National instruments board and in house software. A pattern of nine square 500 micron features was patterned with 20 micron spacing at 99% scan rate with 100 microwatts of power input into the microscope ~ 25 microwatts output from objective.

Unpolymerized monomer was removed by washing with diethylether. The chamber was then filled with a 10 % 1,4-Bis(3-aminopropoxy)butane solution in dimethylformamide (DMF) for 15 minutes at room temperature. The chamber was rinse with DMF and then a solution of 14mg NVOC, 10 uL diisopropylethylamine (DIPEA), and 500 uL DMF was added and allowed to react for 35 minutes. The system was again rinsed with DMF and filled with dioxane. Two squares were scanned on the same laser system described above, with 10 scan lines per feature at 10% scan rate, one was scanned with 5 mW input power and the other 1 mW input power.

The chamber was then filled with a solution of 10mg dansyl chloride, 10 uL DIPEA, and 500 uL DMF and allowed to react for 15 minutes. The chamber was then rinsed with DMF to remove the excess dye and imaged on the same scanning laser system, where the 40x objective had been replaced by a 10x 0.30 NA objective. Fluorescence was collected via an avalanche photodiode which was processed by a Becker and Hickl Time-Related single photon counting module.



The data shows that the feature patterned with 5 mW (white) is significantly brighter than the other features (black/gray):

D. Fill in the following dates:

1. Conception Early March
2. First disclosure to another In early to mid-march I told a lab tech and undergraduate who work with/for me of the concept and told them that I was going to be focused on it from that point forward
3. First written record 03/11/2004
4. First experiment demonstrating the invention 04/14/2004

E. This invention can be used as Enhanced microarray technology where the signal is orders of magnitude larger making it easy to detect binding events (less sensitive instrumentation). Immediate application of this would be for DNA microarrays. This could also be used for the synthesis of large arrays of heteropolymers that could be used in drug development, molecular evolution, or sensor development. Hence, Lab on a chip analytical applications, analytical devices, or microsensor applications. Using multiphoton detection three dimensional surfaces can be decorated with functional groups using this method, these could be cell recognition factors, allowing the construction of complex three dimensional cellular arrays which could be used as bioreactors or artificial organs. These three dimensional surfaces could also be functionalized to form novel biomaterials or drug delivery systems.

F. The problem which this invention solves is Low signal from microarrays which very sensitive equipment to detect. inability to construct three dimensionally dimensionally functionalized materials with a high level of spatial control. Difficulty in screening and or encoding combinatorial libraries.

G. The closest prior art is Frechet Jean M.J. et al, Journal of Polymer Science Part A, 2002, Vol 40, 755-769 and Macromolecules 2003, 36, 1677-1684, used photolithography to prepare monolithic polymers in a spatially defined manner in glass capillaries.

Fodor et al, Science, vol 251, 767-773. Used photolithography in combination with a nitroveratryloxycarbonyl (NVOC) photolabile protective group to synthesize arrays of peptides on a glass substrate.

Satoshi Kawata et al, Nature 2001 vol 412 page 697-698. has created submicron objects using photopolymers in conjunction with two photon excitation.

Shoji Maruo et al, Sensors and Actuators A, vol 100, 70-76. Has used single photon excitation to create 430nm photopolymer features.

H. This invention differs from the closest prior art in that

Frechet used a very different chemistry. He uses free radical polymerization developed by (Ranby, B.) to graft polymers onto polymer surfaces. Therefore he can 'grow' a polymer on the surface of another polymer in a two dimensionally defined way. However, he does not have precise control on the products and can not synthesize heteropolymers with defined sequence in a spatially defined way. Where we propose spatially defined step wise synthesis; protect, deprotect, couple protected monomer/polymer, repeat until a complex structure is created. His method does not have the possibility of doing three dimensional controlled surface functionalization.

Fodor uses glass rather than a patterned polymer surfaces, he does use the photolabile protective in a repeated cycle of coupling and deprotection steps. He is limited to two dimensional patterning.

Kawata and Maruo have constructed three dimensional polymer devices using scanning microscopes but have not, to our knowledge done and synthesis or functionalization on the polymer surfaces.

I. This invention provides the following advantages:

The ability to construct spatially defined functionalized polymer structures of great diversity in two and three dimensions. The larger surface area of the polymer vs. glass increases the signal and potentially the sensitivity vs. conventional microarrays.

Patent Disclosure

04/16/2004

Inventors:

Trent Russell Northen, Neal Walter Woodbury, Sudhir Gudala.
Department of Chemistry and Biochemistry, Arizona State University

Date of invention: 03/11/2004

Date of proof of principle: 04/14/2004

Title: Light directed solid phase synthesis on patterned photopolymers.

Overview:

A new technology has been developed by the inventors that has potential to offer significant advantages over existing microarray technology with long term applications to sensor development, drug development, drug delivery, molecular evolution, and biomaterials.

This invention hinges on the combination of three existing technologies: 1. photopolymers (photoresist – plastic materials that can be patterned in three dimensions with light), 2. photolabile protective groups (these enable the patterned formation of specific chemical bonds in three dimensional space), and 3. solid phase synthesis (the process of generating complex heteropolymers such as, but not limited to, DNA and protein with known sequences on solid surfaces in a completely automated fashion). The combination of these three existing technologies allows for the construction of three dimensional arrays and devices that have tailored chemical functionality. Because of the three dimensional aspect of the photopatterning, an increase of several orders of magnitude in signal strength from, for example, DNA can be obtained as well as dramatic increases in the array density (both two and three dimensional arrays). It has the potential to change the paradigm for the current technology of both DNA and peptide arrays. In the long term this technology could be used to make sensors, combinatorial chemistry platforms, drug delivery devices, biomaterials, and even serve as the scaffold for generating artificial organs.

Huge Opportunity:

"The total biochip market size in 2001 is about \$740 million and may more than triple in revenues, to about \$2.47 billion in 2006" (Technology, Strategic Alliance, Patent Dispute and Market Update – 2002).

Briefly, light is used to make little arrays of solid phase synthesis polymer (which is basically the same thing as photoresist used the electronics industry) on a glass substrate. This polymer can be made in such a way that it is either porous or has a very rough surface (very large surface area) that is covered with reactive chemical groups (for example, primary amines). The reactive groups are then made unreactive by adding a

special blocking group that is photolabile (can be removed by exposure to light). Now, individual polymer elements of the array can be illuminated making them reactive in a patterned fashion and then reacted with materials of interest. One can then add specific molecules only to the elements that have been illuminated. If the molecules added themselves have reactive groups that are blocked by a photolabile blocking group, the process can be continued in layers, building up specific heteropolymers in a patterned fashion. The photolabile blocking group chemistry is the same as what has been used by Affymetrix and other companies to make DNA arrays. The difference is that instead of a monolayer of DNA (or peptide or other heteropolymer) on a surface, one has a much larger number of molecules in the same 2-dimensional element because of the third dimension afforded by basing the array on porous or rough-surfaced polymer elements. This greatly amplifies the signal, making it much easier to detect (the fluorescence from dye reacted directly with the polymer elements is easy to see by eye).

Time Imperative:

There is a group at Berkeley that has recently (2003) published very relevant work. That work does not yet include combining patterning of chemicals on polymer elements using photolabile blocking groups, but they have all of the technology available to them if they were to decide to go in this direction. It is critical that we move quickly to secure our rights to this potentially very significant invention.

Patent discloser:

Inventors: Trent Northen, Sudhir Gudala & Neal Woodbury

Date: 04.16.2004

Title: Light directed solid phase synthesis on patterned photopolymers.

Summary: A general method has been developed to create polymer features and modifying the functionality of the polymer in a spatially resolved manner using photolabile protecting groups to control the addition of desired functional groups.

This invention combines several existing technologies in a novel and useful way. The relevant technologies include: Solid phase synthesis, light directed polymerization, and light directed polymer synthesis.

Disclosed is the method of making polymer structures that have spatially defined chemical features through 1.) photopolymerization to form polymer structures 2.) protection of functional polymer features with photolabile protective group(s) 3.) photodeprotection of desired polymer features, 4.) reaction of deprotected reactive sites with desired reactive chemical species, and 5.) if desired repetition of these steps to form complex functional features.

The features of the polymer and photodeprotected region can be controlled through the modulation of the irradiating light. Small features (~1 micrometer) are generated using high numerical aperture objective lenses and even smaller features can be made using multiphase excitation (50-1000nm) or classical masking methods used in the semiconductor industry.

Background:

Photopolymers:

Photopolymer photo resists are well known and have been used for many years to create small features in the microelectronics industry. More recently they have been used in rapid prototyping or stereo lithography:

- Jan F. Rabek Mechanisms of photophysical processes and photochemical reactions in polymers 1987 John Wiley and Sons Ltd.

Most recently photopolymers have been used in conjunction with high numerical aperture lenses and multiphoton excitation to create very small three dimensional objects.

- Satoshi Kawata and coworkers, Advanced Materials 2003 vol 15, 2011-2014 has used single and multi photon interferential patterning to generate features as small as 50 nm.
- Satoshi Kawata et al, Nature 2001 vol 412 page 697-698. has created submicron objects using photopolymers in conjunction with two photon excitation.

- Shoji Maruo et al, Sensors and Actuators A vol 100, 70-76. Has used single photon excitation to create 430nm photopolymer features.

Spatially resolved biopolymer synthesis is well known and has been used for years to synthesize DNA arrays on glass substrates:

- Fodor et al, Science, vol 251, 767-773. Used photolithography in combination with a nitroveratryloxycarbonyl (NVOC) photolabile protective group to synthesize arrays of peptides on a glass substrate.
- McGall et al, JACS, 1997 vol. 119 page 5081-5090. Used photolithography in combination with the 5'-((α -methyl-2-nitropiperonyloxy)carbonyl) (MeNPOC) to synthesize DNA arrays on glass substrates.
- Michael R. Sussman and co workers, Nature Biotechnology, vol 117, 974-978 used micromirror arrays in conjunction with the MeNPOC protective group to synthesize DNA microarrays.
- Gerard Cagney and coworkers, Nature Biotechnology, vol 18, 2000, 393-397 discusses different applications of protein and peptide arrays.

Solid Phase Synthesis (SPS) is well known and is a method of choice for synthesizing biopolymers (peptides, DNA, etc):

- Merrifield R.B., JACS 1963 Vol 85, 2149-2154 first synthesized a tetrapeptide on a solid resin particle (polystyrene).
- Barany G. et al, JACS 1996, vol 118, 7083-7093 has synthesized a solid phase resin that swells in both water and organic solvents using various methacrylate resins.
- Frechet Jean M.J. et al, Journal of Polymer Science Part A, 2002, Vol 40, 755-769 and Macromolecules 2003, 36, 1677-1684, used photolithography to prepare monolithic polymers in a spatially defined manner in glass capillaries.

Solid phase synthesis techniques have been used to generate combinatorial libraries. These methods have become common to the art, they typically include, dividing the SPS beads into pools after each synthesis step to generate large libraries of peptides. The peptide can be screened and cleaved from the bead can be encoded with some sort of tag for identification

- Lam Kit S. Chem. Reviews 1997, 411-448 this "One-Bead-One-Compound" method.

Photolabile protecting groups:

- Bochet Christain G., Journal of the chemical society, Perkin Transactions 1 2002 vol 2 125-142. Reviews the most common photolabile protective groups.

Biomaterials:

- Langer R. et al, Nature, vol 428 2004 487-492. Reviews biomaterial technology.
- Fisher J.P. et al Annu. Rev. Mater. Res. Vol 31 2001 171-181 describes photoinitiated polymerization and polymer crosslinking for biomaterial synthesis.

One of the significant disadvantages to the existing methods for spatially resolved biopolymer synthesis (Fodor, McGall, and Sussman) is the limited number of reactive sites available on the glass surface (McGall estimates 10-30 picomole/sq-cm). Characterization of reaction products becomes very difficult, requiring sensitive techniques and instruments, for example the most common technique, which is well known to one skilled in the art, for using and characterizing DNA arrays the hybridization of fluorescence probes and use of a scanning epifluorescent microscope to detect these probes. In the case of DNA since a fluorescently labeled complimentary strand can be made for each array element, it would in theory, be possible to characterize any DNA microarray with this technique under the appropriate hybridization conditions.

Since peptides cannot be probed in this same way, due to the non-complimentarity of their structures, other more complicated systems are used. Most commonly, the use of antibody systems in which one antibody is labeled with a fluorescent dye and one antibody (could be the same) is specific for the peptide sequence to be probed (Fodor). This is useful for a proof of principle, but would be impractical for probing large number of peptides.

Even though techniques have evolved to allow the synthesis and screening of libraries using SPS techniques (SPS) screening of the beads is complex.

Polymer structures can be functionalized using this method with cell recognition factors or binding factors such as.....these can be added/deprotected in a spatially controlled manner to create tailored structures. This invention combines the benefits of the array format, large number of reactive sites available in porous solid phase synthesis resin and the ability to form polymer structures using photopolymers. Resulting in larger signals, improved contrast ratios, and better applicability of analytical characterization techniques than existing microarray methods. Since the array is positionally encoded, it is easier to screen and probe than the split pool methods.

These are simply a result of having larger number of sites then on the glass substrate. So that the fluorescence signal is larger when using fluorescent probes or the amount of product produced in a given is large enough to be able to characterize products cleaved off the resin by common analytical techniques such as mass spectroscopy, FTIR, etc.

Further, the array format spatially encodes the peptides so that it is easier to probe than the split pool libraries. These arrays can be probed with analyte for sensor development, drug discovery, or for cell adhesion in biomaterial development.

This invention allows the generation of small three dimensional structures that can be functionalized in spatially defined ways for the construction of sensors, catalysis, biomaterials, drug delivery, molecular evolution, etc.

Summary of the invention:

The system is composed of a photopolymer bearing a reactive group, photolabile protecting group(s), groups to be attached that can also contain the photolabile protective group(s), and devices for illuminating the sample and introducing/removing new reagents. Groups to be attached are not limited to single molecules but could also include macromolecules and even cells.

Polymers/monomers:

For a system where it is desired to detect fluorescence from the array it is important that the polymer system not absorb the excitation light and that it not emit at the detection wavelength. In this case any nonfluorescent nonabsorbing (at the deprotection wavelength) and nonemitting (at detection wavelength) polymer or monomer systems can be used including monomers which are polymerized or polymers that are crosslinked or both. One or more of the following: acrylate, methacrylate, urethane, epoxy, urea, cellulose monomers, protein, glycols, lactic acid, ϵ -caprolactone, trimethylene carbonate, N-vinylpyrrolidinone, 2,2 dimethoxy-2-phenylacetophenone, esters, DNA, RNA,containing side chains.....
or polymers of these monomers and or combinations of these monomers.

Solvents can be incorporated into these systems to modify the pore structure of the polymers. Solvents can include alcohols (methanol, ethanol, butanol, isopropanol, cyclohexanol), acetone, acetonitrile, toluene, etc.

Most preferred are methacrylates and acrylates.

Functionalization:

Polymers/monomers can themselves contain pendent reactive groups like hydroxyls, epoxy, amino, etc groups or they can be incorporated after the polymerization reaction.

Photoinitiators:

Photoinitiators (adapted from JP Fouassier progress in organic coatings vol 47, 2003 16-36) can include in the general classes of initiators: halogens, halogenated organic compounds, hydrogen peroxide, alkyl hydroperoxides, cumene hydroperoxide, peroxides, benzoyl peroxide, non-ketonic peresters, ketones, quinones, polycyclic hydrocarbons, azocompounds, hydrazones, cyclic acetals, 1,3-dithiolane, saccharides, metal oxides, ion pair complexes, metal chlorides, uranium salts, metal carbonyls, metal acetylacetonates, ferrocene, metal complexes, dyes, and polymeric photoinitiators. More specifically radical initiators: azides like azobisisobutyronitrile and derivatives, ketones like benzophenone, thioxanthone, acridone aromatic diketones and derivatives, ketocoumarins and coumarins derivatives; dyes (e.g. xanthene dyes such as eosin (EO) or Rose Bengal (RB), thioxanthene dyes or cyanins); thioxanthenes; bis-acylphosphine oxides; peresters; pyrylium and thiopyrylium salts in the presence of additives such as a perester; cationic dyes containing a borate anion; dyes/bis-imidazole derivatives/thiols; PS/chlorotriazine/additives; metallocene derivatives (such as titanocenes); dyes or ketones/metallocene derivatives/amines; cyanine dyes in the presence of additives;

dyes/bis-imidazoles; miscellaneous systems such as phenoxazones, quinolinones, phthalocyanines, squaraines, squarylium containing azulenes, novel fluorone visible light PIs, benzopyranones, rhodamines, riboflavines, RB peroxybenzoate, PISs with good photosensitivity to the near IR, camphorquinone/peroxides, pyrromethane dye, crystal violet/benzofuranone derivatives, two color sensitive systems, etc.

Colored cationic PIs (such as iron arene salts, novel aromatic sulfonium or iodonium salts) and PS/cationic PI (where PS can be hydrocarbons or ketones or metal complexes) can help to shift the absorption in the visible wavelength range.

Non-ionic photoacids and photobases for the generation of active species in photoresists technology are developed. By now, the design of colored species and proposals of PS for their decomposition remains attractive challenges.

Excited state processes of photosensitive systems for laser beams and/or conventional light sources induced polymerization reactions have been reported in recent works. Typical photosensitive systems under visible lights are classified as One-component system (such as bis-acylphosphine oxides, iron arene salts, peresters, organic borates, titanocenes, iminosulfonates, oxime esters, etc. Two-component system (working, e.g. through electron transfer/proton transfer, energy transfer, photoinduced bond cleavage via electron transfer reaction, electron transfer), Three-component system (where the basic idea is to try to enhance the photosensitivity by a judicious combination of several components).

Most preferred are Azoisobutyronitrile and its derivatives.

Photolabile protecting groups:

Photolabile protecting agents (from Bochet) can include: *o*-Nitrobenzyl alcohol derivatives, α -Ketoester derivatives, Benzophenone reduction, Photosolvolysis-related reactions, Benzyl alcohol derivatives, Benzoin esters, Phenacyl esters, Acylating agents, Fluorene-carboxylates, Arylamines as photo-reductors, Benzophenone as photooxidant, Photoisomerisation *trans*-*cis*, Cinnamyl esters, Vinylsilanes substituted. Most preferred are nitroveratryloxycarbonyl, 5'-((α -methyl-2-nitropiperonyloxy)carbonyl)

Groups to be added:

Groups to be added onto the polymer structures include, sugars, amino acids, nucleic acids, multifunctional amines, ethylene glycol, acid labile groups, base labile groups, dyes,and combinations of or polymers of these monomers. Sequential light directed synthesis can be used to build complex sequence specific polymers.

Most preferred groups include amino and hydroxyl groups.

Method of light modulation:

Light can be modulated (spatially patterned) using a scanning laser system composed of a laser, shutter, microscope objective and stage. In this case the stage movement and

shutter are controlled so that the shutter is only open when the stage is positioned so that the light will illuminate a desired position.

Photolithography is well known to the art but briefly it utilizes masks where light is blocked by some parts of the mask and not others. In this way the illumination reaching the sample can be controlled. Light sources typically include lamps or lasers.

Micromirror arrays are a more recent way of modulating light. By changing the angle of the mirrors in the array light can be directed towards a surface or not. In this way light from an excitation source (lamp or laser) can be selectively reflected onto desired regions of the sample to be exposed.

The preferred embodiment is either a micromirror array or scanning laser system

Substrate: Substrates can include glass, quartz, silicon oxide or other metal oxide surfaces, polymers bearing reactive groups. It is not necessary that they be transparent since illumination can be from above. In the case of glass, quartz, and silicon oxide these surfaces can be modified to react with the polymer for a covalent linkage, though this may not be desirable or necessary in all cases since intermolecular attractive forces can be used to 'glue' the features to the substrate. Where modification is desirable silanes common to the art can be used, the most common being aminopropyl triethoxysilane or 3-(trimethoxysilyl)propyl methacrylate.

The preferred embodiment is glass cleaned with acid and base as described in McGall JACS 1997 and functionalized from a 1% solution of 3-(trimethoxysilyl)propyl methacrylate in 95% ethanol 5% water.

System for introducing reagents: Systems for introducing and removing reagents include an optical flow cell coupled with manual or automated introduction and removal of reagents. Wells or plates where reagents are introduced manually or by automation. Automation is provided by machines such as peptide synthesizers that are designed to introduce and remove reagents.

Analytical Techniques:

Array elements can be probed in situ through various spectroscopic techniques including fluorescence, absorption, infrared spectroscopy, raman spectroscopy, nonlinear spectroscopy, and surface plasmon resonance or elements can be removed from the surface through the use of labile linkages between the coupled material and the polymer. Thus the material can be cleaved and a host of analytical techniques can be used including HPLC, NMR, Mass spectrometry, capillary electrophoresis.

Most preferred include fluorescence detection of hybridized, bound, or covalently linked probes or groups, infrared spectroscopy, and mass spectrometry of cleaved materials.

Example 1:

A glass cover slide was cleaned for 15 min at RT with 60/40 sulfuric acid/hydrogen peroxide, placed in 10% sodium hydroxide at 70 C for 3 min, placed in 1% HCl at RT for 1 min, between each step it was soaked in nanopure water for 3 minutes. A solution of 1% 3-(trimethoxysilyl)propyl methacrylate in 95% ethanol 5 % water was made and mixed for 10 minutes, the slide was then added and left to react at RT for 15 minutes with gentle agitation. This slide was soaked in isopropyl alcohol for 3 min then nanopure water for 1 min then placed in a 100 C oven for 5 minutes after which the oven was turned off and nitrogen was blown through for 1 hr. The slide was stored under nitrogen until it was used.

A blend of methacrylate monomers and photoinitiator (900 uL trimethylolpropane trimethacrylate, 100uL glycidyl methacrylate, 10 mg azobisisobutyronitrile) was prepared and nitrogen was bubbled through for 15 minutes before it was injected into a Focht cell (Biopetech inc. Butler, PA) that had been flushed with argon. The above slide was mounted in the flow cell. This cell was then mounted onto a Prior scientific microscope stage on a Nikon microscope and illuminated through a 40x 0.75 NA objective with 370nm light from a Ti:Sapphire laser modulated by a Conoptics modulator (shutter). The system was controlled via National instruments board and in house software. A pattern of nine square 500 micron features was patterned with 20 micron resolution at 99% scan rate with 100 microwatts of power input into the microscope ~ 25 microwatts output from objective.

Unpolymerized monomer was removed by washing with diethylether. The chamber was then filled with a 10 % 1,4-Bis(3-aminopropoxy)butane solution in dimethylformamide (DMF) for 15 minutes at room temperature. The chamber was rinse with DMF and then a solution of 14mg NVOC, 10 uL diisopropylethylamine (DIPEA), and 500 uL DMF was added and allowed to react for 35 minutes. The system was again rinsed with DMF and filled with dioxane. Two squares were scanned on the same laser system described above, with 10 scan lines per feature at 10% scan rate, one was scanned with 5 mW input power and the other 1 mW input power.

The chamber was then filled with a solution of 10mg dansyl chloride, 10 uL DIPEA, and 500 uL DMF and allowed to react for 15 minutes. The chamber was then rinsed with DMF to remove the excess dye and imaged on the same scanning laser system, where the 40x objective had been replaced by a 10x 0.30 NA objective. Fluorescence was collected via an avalanche photodiode which was processed by a Becker and Hickl Time-Correlated single photon counting module.

The data shows that the feature patterned with 5 mW is significantly brighter than the other features (see image below):

Example #2

A monomer mixture of the following composition was prepared: 1 mL hydroxyethyl methacrylate 2.6 mL Trimethylol propane trimethacrylate and 36 mg azobisisobutyronitrile. They were then sonicated for 5 minutes. Nitrogen was bubbled through the sample for 5 minutes. A biopetecs FSC2 chamber was purged with argon and then filled with the nitrogen flushed monomer mixture, with coverslip functionalized with trimethoxysilyl propyl methacrylate as described in example #2. The chamber was mounted on modified Prior scientific Proscan stage attached to a Nikon microscope. Laser excitation was obtained from a Spectra-Physics Tsunami mode-locked Ti:sapphire laser (742nm), which went through a Conoptics shutter and was later doubled to 371nm. This then was focused through a Nikon 0.30 NA 10x objective onto the sample. Laser power was set to 250 microwatts going into the laser with approximately half of that power at the sample. The photopolymer was patterned using in-house software, designed to control the stage and shutter. The following patterns were made:

Position	Power into laser (μW)	Exposure time (ms)	Number of features	Feature spacing (μm)	Focus vs. cover glass (μm)
1	500	500	5x5	1000	250 below
2	250	1000	5x5	1000	250 below
3	250	1000	5x5	1000	250 above
4	250	250	10x10	250	250 above
5	250	250	10x10	500	250 above
6	250	250	20x20	250	250 above

After patterning the chamber was rinsed several times with dimethyl formamide (DMF) to remove unpolymerized monomer. The hydroxyl groups of the polymer were coupled to a NVOC protected glycine (NVOC-gly). This was prepared from the amino acid and NVOC-acid chloride using schotten bauman procedure. The NVOC-Gly was activated with the coupling agent O-Benzotriazole-N,N,N',N'-tetramethyl-uronium-hexafluorophosphate (HBTU) using the following procedure: 21 mg NVOC, 24 mg HBTU and 1 mL DMF were mixed and allowed to react for 30 seconds and 12.4 μL Diisopropylethylamine (DIPEA) was added, this mixture was allowed to react for 3 min before adding to the chamber. This mixture was allowed to react in the chamber for 30 minutes without mixing and then another 30min with recirculation.

Any unreacted sites were acylated with acetic anhydride: a solution of 5 mL DMF, 146 μL DIPEA, 100 μL acetic anhydride was prepared. The chamber was flushed with DMF and then filled with this solution and allowed to react for 30 minutes. The chamber was then flushed with DMF and then filled with dioxane.

Features were deprotected using the same apparatus used to pattern the photopolymer. The laser was set to the same wavelength and the power was set to 500μW. For the larger patterns (1, 3, and 4) by manually finding the feature to be deprotected using the microscope and opening the shutter for ~ 30 seconds to expose the feature. The smaller

patterns (5 and 6) were deprotected by scanning adjacent features with the laser beam in a series of parallel lines forming squares. Position #5 was scanned with six 700 μm squares with 700 μm spacing between squares and each square was composed of 30 scan lines scanned at 1% of the maximum scan rate and 1 mW power into the microscope. Position #6 was scanned with 312 μm squares with 448 μm spacing with 20 scan lines at 1% of the maximum scan rate and with 1 mW power into the microscope.

After scanning the dioxane was drained and the chamber was flushed with DMF. The chamber was then filled with 10 mg/mL fluorescein isothiocyanate (FITC). This was allowed to react for 15 minutes. The chamber was flushed with DMF and sat overnight filled with DMF. The chamber was then flushed with fresh DMF and imaged.

Imaging was performed using the same scanning laser apparatus and laser configuration as used for the patterning and deprotection of the polymer-NVOC-GLY. The laser was set to 8 μW input power into the microscope. Emission from the FITC was collected by an APD detector. Imaging was done with in-house software using a Becker & Hickl GmbH SPC-830 high performance photon counting board. All images were 64x64 pixels with the stage moving at 10% of maximum scan rate.




Image 1: Position 1 Shows the ordered spacing of the features (dark spots are unpatterned bright spots and lines have FITC. Note that the polymer is long hair like structures, some of which have fallen over.




Image 2: Intensity x cross-section at scan pixel 23 of position 1. Note that the dip at $y=24$ corresponds to a unpatterned feature and the peak at $y = 46$ corresponds to a patterned feature. The image scan spacing was $50\text{ }\mu\text{m}$ so the two features are $\sim 1\text{ mm}$ apart which corresponds to the distance between features in pattern 1. The contrast ratio is very high since the unpatterned feature is darker than the background.




Image 3: Fluorescence intensity from position 3 shows alternating features of protected (dark) and deprotected (yellow, green, pink).




Image 4: Intensity x cross-section at scan pixel 30 of position 3. Note that the dip at $x=15$ corresponds to a unpatterned feature and the peak at $x=2$ corresponds to a patterned feature. The image scan spacing was $75\text{ }\mu\text{m}$ so the two features are $\sim 1\text{mm}$ apart which corresponds to the distance between features in pattern 1. The contrast ratio is very high since the unpatterned feature is darker than the background.




Image 5: Fluorescence intensity from position 4 showing the alternation of protected and deprotected features.




Image 6: Fluorescence intensity from position 5 showing the box like deprotection pattern.




Image 7: Fluorescence intensity from position 6 showing the box like deprotection pattern.





Image 8: Fluorescence intensity from position 6 showing the box like deprotection pattern of small features note that the features are 4 pixels in diameter and 25 pixel spacing, with a 10 μM scan spacing that is 40 μM and 250 μM spacing as expected.



Claims:

1. A general method to create spatially defined complex polymer structures through sequential deprotection and addition of polymers/monomers to photopolymer structures.
2. A way of generating polymer arrays using photolabile groups with acrylate and methacrylate monomers and AIBN and its derivatives as a photoinitiator.
3. A method of enhancing the signal from microarrays by constructing the microarrays on patterned photopolymer arrays.
4. A method of generating three dimensional structures functionalized with a plurality of spatially defined functional groups.

Materials:

02/20/2004

RenShape SL 5510 SLA System From Vattico

Propylene Carbonate laser 10mw or more @ 100 rpm

Isopropyl Alcohol laser my w/ N₂

T: Sapphire laser 10 in FWHM 80mHz 730nm / 2

10x 0.7NA objective


Experiment	Scan Speed	Line Spacing	Power	Size	Obs
1 'ASU'	1%	40um	1uW	5.5 x 5.5mm	1/2
2	1, 3%, 100%	40um	100uW	5.5 x 5.5mm	1/2
3	1, 2, 3, 4%	50um	100uW	4 @ 1um	1
4 'ASU'	1%	40um	100uW	6 x 6mm	Take a long time
5	Repeated #4				
6 (2nd 'ASU')	1%	9um	10uW	500um x 500um	
16	500 x 500um	1%	50um	10uW	(4500um) ²
16	500 x 100um	1%	10um	10uW	(900um) ²

* IDEA ① Use polyacrylonitrile attach pellets
 ② or ~~epoxy~~ epoxy / methacrylate & decure? attach pellets

make an array of PADS laser gate of pulsed
 laser sensors

Continued on Page

Read and Understood By



Signed

4/16/04

Date



Signed


4/21/04

Date

Purpose: #1 to determine the threshold
for photo polymerization of 2nd
excitation @ 730nm

→ 11nm FWHM

Laser 730nm output after Beam splitter before shutter = 650
Output from objective 280mW 0.8mW



'Grayscale M'

400nm x 400nm

8 bit RLE


X step = 2nm

Y step = 6nm

Scan rate 2638 51 63 76

This worked but the pattern came off when I
washed / wiped it may have been focused
above the surface.

Laser
output
before
shutter
2nd



Laser at 240mW

190nm x 650nm

10nm x steps

10nm y steps

10x objective

US22 Pattern Response 1/13/04 as a substrate

There was some polymerization

only near where there is significant polymerization

is between areas (where the laser is changing scan speeds)

Need to systematically slow scan speed

There was a visible border where I lost more light

Continued on Page _____

Read and Understood By

Signed

4/16/04

Date

Craig L. Mager

Signed

4/21/04

Date

Process

See if it is possible to couple a Diels-Alder to a reactive epoxy resin on the 5510 polymer

213 mL 83% pure Propoxy Butene + 30 mL DMF

Scanned $6 \times 9 (250 \text{ nm})^2$ films w/ 200 nm spacing 19 mm

Spacing @ 2 mm w/ 40x objective 50% SS

Washed w/ propylene carbonate for ~ 3 min

Washed w/ 10% NaOH 2. Wash w/ N_2

Washed to Diels-Alder solution @ 4:15 PM

Scanned same but w/ 100 nm intensity @ 100% SS


Still looks like too much

Try 10 mW for 1st 30 sec

Switched to 1 mW for the next

threshold

@ 10 mW
Made reactions
like this



didn't see anything @ 1 mW

need to either start with polymer or use two photos
overcoat

- Sample left in Diels-Alder solution at glass
residue a small piece & washed w/ methanol
& placed in light - very soft until cured

- may want to cure while reacting w/ Diels-Alder

6:55 PM - Added two samples (1 w/ & 1 w/o Diels-Alder) to 100% (w/ Diels-Alder)
8:00 PM - Added 60 mL DMF

Washed 50 mL methanol

+ Diels-Alder is unbelievably fluorescent - not very, if at all fluorescent

so this may be a possibility

epoxy has good resistance to
solvents (hexane,
DCM, THF, Acetone, etc.)
Continued on Page

Read and Understood By

Not DME

Signed

Date

Signed

Date

PROJECT SAW in a lab at 400

Notebook No. _____
Continued From Page _____

Start w/ SAW input ~ 1mW output from 400. 1J/cm² @ 88% SR
Start High & Go Low Curves measure of black growth
1st Scan 21

<u>0.25 J/cm²</u>	<u>0.125 J/cm²</u>	<u>0.0625 J/cm²</u>	<u>0.03125 J/cm²</u>
#5	Focused w/ 4 of 4.6	Scanned P-35. White coat!	#3
4000, 0	4000, 3000	4000, 6000	4000, 9000

<u>9 J/cm²</u>	<u>3.5 J/cm²</u>	<u>1 J/cm²</u>	<u>0.5 J/cm²</u>
#1	#2	#3	#4
0, 0	0, 3000	0, 6000	0, 9000

1st Scan	300x300 w/ 150 spacing	20 Scan lines	3x3	SR=10%	9 J/cm ²
2nd	Same			SR=25%	3.5 J/cm ²
3rd	Same			90%	1 J/cm ²
4th					0.5 J/cm ²
5	300x300	Same		25mW input	90% 14mW
6				1.25mW input	90% 0.25 J/cm ²
7				625mW	90% 0.125 J/cm ²
8	200x250	50 lines	5mW spacing	313mW	90% 0.0625
9	Same			15L	0.032

* More of Sigmahurt Scatter off top low 51/!

Added dye to y/m input + saw spec but not for 20 min

* Focusing is very difficult
Went up ~ 100-180 microns
By 45X on the size
[Signature]
Signed _____ 4/16/04 Date

Read and Understood By
Wish felt it was the low
Craig Chayer
Signed _____ 4/21/04 Date

Continued on Page _____

2
250
550
441
100
90

Inching of 545 RF 35

w/ 50000 Power input

Dist is too pleasant

Inches (200-1600)² 15x15lots of present (fairly sharp Co₂H₂) on surface

Inches or more

x=12

Between future & 2000

Dist. from

Inch intensity

Background intensity

Comments

Inches (S)

100-550

100-250

100

90

110

120

130

140

150

160

170

180

190

200

210

220

230

240

250

260

270

280

290

300

310

320

330

340

350

360

370

380

390

400

410

420

430

440

450

460

470

480

490

500

510

520

530

540

550

560

570

580

590

600

610

620

630

640

650

660

670

680

690

700

710

720

730

740

750

760

770

780

790

800

810

820

830

840

850

860

870

880

890

900

910

920

930

940

950

960

970

980

990

1000

1010

1020

1030

1040

1050

1060

1070

1080

1090

1100

1110

1120

1130

1140

1150

1160

1170

1180

1190

1200

1210

1220

1230

1240

1250

1260

1270

1280

1290

1300

1310

1320

1330

1340

1350

1360

1370

1380

1390

1400

1410

1420

1430

1440

1450

1460

1470

1480

1490

1500

1510

1520

1530

1540

1550

1560

1570

1580

1590

1600

1610

1620

1630

1640

1650

1660

1670

1680

1690

1700

1710

1720

1730

1740

1750

1760

1770

1780

1790

1800

1810

1820

1830

1840

1850

1860

1870

1880

1890

1900

1910

1920

1930

1940

1950

1960

1970

1980

1990

2000

2010

2020

2030

2040

2050

2060

2070

2080

2090

2100

2110

2120

2130

2140

2150

2160

2170

2180

2190

2200

2210

2220

2230

2240

2250

2260

2270

2280

2290

2300

2310

2320

2330

2340

2350

2360

2370

2380

2390

2400

2410

2420

2430

2440

2450

2460

2470

2480

2490

2500

2510

2520

2530

2540

2550

2560

2570

2580

2590

2600

2610

2620

2630

2640

2650

2660

2670

2680

2690

2700

2710

2720

2730

2740

2750

2760

2770

2780

2790

2800

2810

2820

2830

2840

2850

2860

2870

2880

2890

2900

2910

2920

2930

2940

2950

2960

2970

2980

2990

3000

3010

3020

3030

3040

3050

3060

3070

3080

3090

3100

3110

3120

3130

3140

3150

3160

3170

3180

3190

3200

3210

3220

3230

3240

3250

3260

3270

3280

3290

3300

3310

3320

3330

3340

3350

3360

3370

3380

3390

3400

3410

3420

3430

3440

3450

3460

3470

3480

3490

3500

3510

3520

try using more concentrated ammonia in TAC. For a short

Feb-09 426 μ L Bis-Ammonia Propylamine in 1 mL TAC
 ~ 180 mM

first scanning at 1.0 mW	$10 \times (2000 \mu\text{m})^2$	25×25	works ✓
1 ps exposure			
scanning	0.5 mW	same	✓
	0.25 mW	same	✓
	0.1 mW	same	no

the stresses were too high! they fell atop lechite

Tried at $40 \times (2000 \mu\text{m})^2$ 20×20 lines

they don't stick to the glass very well

may also be APPS surface


keep working at $10 \times$... can always (hopefully) scale down

250 μ L $10 \times$ 1 ms / pixel $(4000 \mu\text{m})^2$ 20×20 pixels
 200 nm between spots

put on PC Agitated gently until excess had dissolved
put microscope lamp on it for 5 mins while in the
PC

Altogether seems to remove 5510 very well
 125 mW at 90% SR $10 \times$ is light at the threshold
only the edges polymerized.

gives us more intensity



Tried making some little 'Hem-1'

Continued on Page _____

1st were @ 0.2 mW Read and Understood By

2nd were @ 0.3 mW

Signed _____

Date _____

$1000 \mu\text{m} \times 4000 \mu\text{m}$ 25×25 lines

Worked with 800 μ m
through 75
Signed _____
4/21/09

PROJECT 3/15/04

Notebook No. _____ 73
Continued From Page 71

Making 'waves' from pg 71

w/ DAWG filters exactly @ 365nm Get a 'huge' fluorescent signal. Probably from photo initiators.

Need to take an excitation / emission spectrum from polymer to figure out if I can find a dye to use with it / where it will have a small background. Also see that it is generating white light.

Making w/ 40x 40xw features are 25um in diameter. Have a S/N of 80 Much history 1160 = 5 N = 38 ^{Dark Count = 5}

What would happen if you scattered blue into there. You would have a SPR system?

$$\frac{1160}{13} = 9$$

Power dependence of emission


Input Power (mw)	APC	Input Power (mw)	APC
0	7.5 E 2	90	7.6 E 6
2	7.8 E 2	80	2.5 E 6
3.4	8.7 E 2	70	2.3 E 6
5.4	9.8 E 2	60	2.2 E 6
8.6	1.1 E 3	50	2 E 6
11.2	1.2 E 3	40	1.8 E 6
16.6	1.4 E 3	30	1.5 E 6
23.4	1.7 E 3	20	1.1 E 6
39	2.3 E 3	10	6.4 E 5
		9.7	5.7 E 5
		8.3	4.8 E 5
		3.6	2.9 E 5
		1.2	1.1 E 5
		.8	7.2 E 4

Not well coupled

Continued on Page

Read and Understood By

1.5 E 3



Signed

4/16/04

Date

Craig Emmer

Signed

4/21/04

Date

3/16/04

Notebook No. 30-211
Continued From Page 1

Continued From Page

[illegible][illegible]

PROJECT

Notebook No.

Continued From Page

Exch. @ 550	Exch. @ 500	
450 2423	450 3341	
475 2966	475 5039	
500 3567	430 2470	
425 2100	420 3055	
400 1781	410 5234	
365 1080	400 3437	500
350 515	365 2771	420
	350 1282	460

100 mm 8.45 Shift

Exch.	Exch. @ 500	Exch. @ 500
350 450	1275	
400 478	2913	
425 500	2278	
450 527	2177	
475 540	1492	
500 534	2330	
365 447	2726	


500

nothing just scatter

(515) may be the music execution number. Need Argon
low level

00 = 3 @ 375 mm
off scale @ 365 mm

Argon line



Continued on Page

Read and Understood By

4/16/04

Craig Mager

4/21/04

Signed

Date

Signed

Date

PROJECT Alkyne fibers out of 5510 w/ 40x

Notebook No. _____


Continued From Page _____

13 μ W input w/ 40x objective 36.5 μ m resolution


APTS SLIDE

3000 x 3000 μ m30 μ m x 30 μ m

Shutter open 1 min



using one to estimate



Proceed into bed of acetonitrile

need for to see thru put over a light (acetonitrile)
for 5 minFibers are in 5 μ m wide 2 μ m long.

Try to make longer fibers by using thicker
 coating placed 2 150 μ m (can slip on APTS slide
 & used a 3" x 2" CS to spread ~ 40% of polymer on
 surface

30 μ W power (1000 μ W/cm²) 25-28 μ m is long / pixel40 μ m spacing

Watched well except fibers fell over ~ 80 μ m long 1.25 μ m
 Continued on Page 77

Read and Understood By.

Signed

4/16/04

Date

Signed

4/21/04

Date

start - SP by coloring in place stay low intensity
for a long time.

137/min for 2hr w/ 500-500

From page 71 rather than high intensity

300 mW for 1 ms try (lower for 30 ms)

if this doesn't work try

flashing twice; once to pattern & the beam time to
cure

☐ 3 ms

☐ 300 ms

Found using mes laser.

☐ 30 ms

10 mW 10x objective $(1000 \text{ nm})^2 / (10 \text{ pixels})^2$

30 ms @ 300 mW, 3 ms

3. (only one that worked, thought it still works
but starts up.

Cuts

Features were ~ 25 μm in diameter

masked w/ 100 mW for 1 ms w/ 40x

worked well!

2 sec @ 10 mW

$(1000 \text{ nm})^2 / 10^2$

} But appeared to

3 ms @ 100 mW

3x

$(800 \text{ nm})^2 / 8^2$

} work

Don't start up as well.

Masked w/ 40x: 40 mW too much power

Very large fibers ~ 2200 μm long 20 μm w.x

6. make 400 x 400 μm / 30 x 30 @ 10 mW w/ 40x

~ 20 μm across ~ 200 μm long

Continued on Page _____

Read and Understood By

This is how to how done the

Bum 13x 4/16/04

Craig Enage

4/21/04

Signed

probably 10 mW & 150 μm is how
a couple.

Date

Signed

Date

New try w/ 40x

$$\frac{\pi}{4} (12 \times 10^{-4} \text{ cm})^2 \quad \frac{18}{5} \text{ J/cm}^2 \quad \frac{26.5}{5} \text{ J/cm}^2$$

was 800ms w/ 2.1mw input w/ 40x
 send so w/ 8.2mw input

The nics method of blowing works well: looking for
 the maxima closest to the surface
 everything has collapsed!

STANARD seems to work well for removing polymer


30W	40x	1ms	200ms	0	0
40	0	1ms	30ms	1000	500 x 500 nm / 10x10 pins
		1ms	100ms	2000	0
		10ms	10ms	3000	0
		100ms	100ms	4000	8000
		100ms	0	2000	2000
		200ms	0	3000	8000

Send the 5510 by blowing Argon across the surface.

gathering this

only the low! low! low! bracket?

Still feel our interest by our pattern come out
 as turning



Continued on Page _____


Read and Understood By

Signed

Date

Signed

Date

P26₁₀ MA

EGDM

1 st Attempt	1 st Attempt	2 nd	3 rd	4 th	5 th
P26 ₁₀ MA	0.9 ml	450 ml	45 ml		
EGDM	19 ml	10 ml	10 ml	100 ml	100 ml
AzoIBN	11 mg	0	2		10 mg
Ethanol	1 ml	0	6		50 ml
Benzophenone	0	6 mg	28 mg	50 mg	
	Mixing	Mixing	Some later formation?		worked after a while

NEED to think about kinetics... NEED

TO REMOVE INHIBITOR!

NEED Anneal Triethylamine @ 50 mM

1st Initiation2nd3rd

50% monomer

	Unit	1	2	3	4	5	6	7	8
EGMA	50-98%	98	98	98	98				
AzoIBN	1%	1	0	0	0				
PB	1%	0	1	1	1				
TEA	5-12	0	0	1	0				
BAPB	5-12	0	0	0	1				
ETO4	0-48%	0	0	0	0				

Continued on Page

Read and Understood By

Signed

Date

Signed

Date

removed inhibitor by running over Sigma-Aldrich
column.

first washed column w/ 2 to 4
vol of solvent & discarded 1st part

flushed w/ N_2

kept under N_2
polymerizing in beaker

$$\left(\frac{50 \text{ mmol BP}}{L} \right) \left(\frac{182.2 \text{ mg}}{1.22 \text{ mmol}} \right) (0.001 L) = 9.11 \text{ mg} \sim 1\%$$

① 10 mg Azobisisobutyronitrile + 100 μ l cytolene grade DMF + 900 μ l MMA

② 10 mg BIP + 10 μ l TCA + 1 μ l MMA

200 μ l ① + (10 mg polyvinyl pyrrolidone + 100 μ l DMF)

200 μ l ② + "

The samples w/ PVP worked significantly better than
those without.

Made $\frac{200 \text{ mg}}{200 \mu\text{l}}$ Azobisisobutyronitrile in DMF $\frac{200 \text{ mg}}{200 \mu\text{l}}$

Made $\left(\frac{10 \text{ mg PVP}}{200 \mu\text{l MMA}} \right)$


added 50 μ l of Azobisisobutyronitrile solution

worked ok, not sure if it is any better than 1%
Need to check even more.

Need controls & other stuff what is happening

Continued on Page _____

Read and Understood By



Signed

4/16/04

Date

Craig Chaga

Signed

4/21/04

Date

	MAR	DMB	PSP	BL	A20	TEA
1		X	X			
2		X	X		X	
3	X					
4	X				X	
5	X		X		X	
6	X	X		X		

#5 10mg PSP + 10ml of 20

A20 + 0ml + 180ml MMT

Exposure 2nd input for 3 min w/ 10x

OBSERVATIONS:

Before wash off excess #4 & #5
have cloudy material on top. The others
are clear

4, 5, 6 are the only ones that polymerized
4 is white!
5 is clear large lumps very nice
6 is very small & clear.

Formula Number 5 seems to be a good
starting point.

Exposure test with Formula #5 2.1 min w/ 10x

I see Bumpy visible spec

10 see Hair < 100 mm wide probable close to 20

30 see ~ 200 mm

40 see 1 mm wide

Good enough for now

Continued on Page

Read and Understood By

Signed

Date

Signed

Date

Try adding diamine 1%

① 10mg BP + 2mg PZDA (1500g/mol) + 10ml DiAzo - DMF + 178μL MMA

② " 2mL 95% ethanol (propoxy Butene) "

Patterned ① : ② : #5 from before 2.1 mW 10x 30 sec

whos w/ acetone


all spots look the same

Stained w/ 10mg DAPI in 500ml DMF + 10ml PZDA
430nm left 1 of two #① w/ stained one on edge
445nm

See figures, but with IAZ #5 + PZG + DAPI 1 of
#5 + DAPI are fluorescent #5 + PZG is
1/10 as fluorescent (excitation @ 365nm, 300mW 400)

thus maybe ok. Need to determine if the dye
is reacting w/ surface or is simply being adsorbed.

Test By Acetylating the surface & staining w/ dye



exposed to 2.2 mW 10x for 3 sec

+ PZG Hx as 1-11 20 - PZG

Ultraviolet light source 515nm: 0.9ml DMF + 100μL AA + 100μL PZDA

Dissolve polymer (cured) w/ 15 min

cured polymer - allowed to dry overnight before dissolve.

& it swells in DMF - adds to remove dye.

needs to cure under lamp for ~ 30 min
to crosslink.

Read and Understood By

Signed

Date

Signed

Date

4/16/04

Craig C. Magill

4/21/04

Continued on Page

Run GMA on Column 18" x 0.5"

1st run 10 ml on F

2nd 10 ml 200X

3rd 10 ml 200X

Con GMA, do not polymerize w/ 1% Azo

Run A Successive Finer (re run)

	% w/m	mm	GMA	PEB DMA	Result
1	X	X			works turned white?
2	X			X	Clear w/tes (30 sec)
3	X	20		40	worked well fast! wh
4	X	30		20	
5	X		80	40	Slow (60 sec) Clear
6	X		X		Nothing
7	X	20	5	75	similar to 9
8	X	30	20	30	very slow (slowest)
9	X	35	5	60	worked the best. skin very slow
10	X	75	5	20	turned white @ 30 sec

Added 10% BAPS in DMF to #7 w/ 10 ml had

formed w/te Gel prod. Most was still w-acted, probably dissolved in DMF & reacted w/ P1Amine.

Come over here for 15 min + report ✓

Reported #9 - 10 min for 5 min, no apparent change in physical properties

Mailed w/ Dwyer (#8 & #9) for 15 min

Continued on Page _____

Read and Understood By _____

Read and Understood By _____

2X w/ Dwyer

Signed _____

Date _____

Signed _____

Date _____

4/21/04


This is good & 9 looks like a good female, may want to
 let her GMA to triple check that the inhibitor is removed.
 There needs a way to attach to the alkyl. & need
 a way to confide

Make more of #9 using GMA that has been
 stored w/ Aldrich beads for 1 hr

Don't like how the chem. turns when
 try w/ MMA

Azo	GMA	per DA
100ml	100ml	800ml
100ml	50ml	850ml

oops ended up adding 2% Azo



1 2 20ml of each.

Has become fluorescent not after
 P26-GMA has become fluorescent! probably from both
 Go see it with stand in. everything mixed together
 messy pattern

Added 10% BAPB 1253 pm

Try 50ml Azo + 100ml GMA + 850ml MMA

Continued on Page

Read and Understood By

Signed

4/16/04
 Date

Craig Magee
 Signed

4/21/04
 Date

→ Purchased more PEGDMA over Sigma tested
 → Recrystallized Azobisisobutyronitrile from MeOH @ RT (AIBN)

Recd ① Some AIBN + 800 μ l PEGDMA + some GMA

② " " + some GMA + 900 μ l PEGDMA + 700 μ l MAA

Some AIBN + 800 μ l PEGDMA + some GMA + some BAPB

→ Cured for 60 sec @ 2.5 mW

→ Measured both w/ DSC both were around #2 more than #1

→ This makes a soft polymer, probably the BAPB is typical when the measured -

→ The MAA seems to make it cross faster & more rigid

→ May be that in very way too much light & getting short polymers.

try less light & blendy w/ some MAA

Use 200 μ l for 60 sec w/
 1% AIBN + 1% (BAPB + GMA) + X

X = 10 - 60% MAA ; 0 - 70% GMA

Continued on Page _____

Read and Understood By

Signed

Date

Signed

Date

Real work of ~~the~~ GMA through Alorich resin

Dispense column (0.5" x 8")

dispensed w/ from resin for PGE DMA

mass 12 AIBol (recrystallized from MeOH) solvent.

blended ^(20ml) N_2 through for 2min

placed on cold slide w/ cold petri dish in and

w/ N_2 blowing through. oven at $85^\circ C$

11:45 am 4/6/04

12:45 pm Both samples cured into clear hard polymers

The PGE doesn't stick to the glass very

well can't cut the GMA off the glass


All these are like little irregular

leaves

laser @ 365nm 2min input into processor

Went sample to 30% over to me I poured w/ N_2 for

twenty minutes laser found to 500-1000nm



1, 2, 3, 4

1 = 30 sec @ 2mw

2 = 30 sec @ 200mw

3 = 300ms @ 200mw

4 = 4 @ 2mw

Continued on Page

Read and Understood By

Signed

4/16/04

Date

Signed

4/20/04

Date

None of the light induced polymerizations today worked very well. The GMA sticks to the Glass like crazy. The PGE doesn't appear to swell very much in DMF or water.

together they might work. Try various concentrations in DMF

PGE GMA

100

0

have been stick to glass

0

60

stick to glass amazingly well

①

50

50

②

99

1

then sparged it w/ for 1 hour

③

90

10

then spalled on plate in DMF @ 20

④

75

25

sparged w/ N_2 

flushed in 15 min.

None stuck

All are clear

At 1st set in BAPB (mt) 62 in 4/6/04

only takes 10 min but need

The 90/10 cracked overnight, may not have anything to do w/ binder

boasted
The set treated w/ Diamine is very fluorescent
even after washing 2x w/ Acetone 30 w/ DMF
The set not treated w/ Diamine is not fluorescent

Continued on Page _____

Read and Understood By

Signed

4/16/04

Date

Craig Magee

Signed

4/26/04

Date

w/ 120 AIRN

	<u>TMA</u>	<u>PEOMA</u>	<u>GMA</u>	
A	100	0	0	TMA 378.4 mg/m ²
B	90	0	10	PEOMA 198 mg/m ²
C	45	45	10	GMA 172 mg/m ²
D	49	0	1 3.5	

A & D are very brittle, they crumble

B & C are much better

B + 10% Air in AME 220 pm

E 75

f

Regressed 4 zone spots of B in area (1A) for 30 min

350w tube for 5 min (300w) The plate

is slide 4 zone spots for 5 min

put in plate dish on plate dish w/ plate

Cured w/ black cloth

Sawed ~~some spots~~

30 min deposed

Set out 5 min deposed

300w 50w
100w 50w
300w 50w
worked

50
50
0

1 sec @ 300w work both (work)

1 sec @ 300w do not work

300w 30 sec (work)

75% SS @ 300w (500x500mm²) 12 min

300w do not work
50w do

appears that you need a lot of light!

Continued on Page _____

Read and Understood By

Signed

04/16/19

Date

Craig L. Magrath

Signed

4/21/04

Date

Take a methacrylate slide put some of the polymer on it
Scanned area at $(1\text{mm})^2$ 1000 lines @ 3 mW
all @ 3.2 mW

○	$(1\text{mm})^2 / (100)^2$	99.2%	SS	@ Focus
○	$(1\text{mm})^2 (50)^2$	99.2%	SS	
○		50%	SS	
○		25%	SS	

Made more B 900 ml TMA + 100 ml GMA + 10 mg AIBN
sargal - 1 ml for 5 min. Gas N_2 through over.
placed slide in over when
poly. done

04/13/04

TMA

2-aminethyl

AIBN

~~Methacrylate~~

HCL

my
MT₁ soluble

50

50

12

90

10

22

95

5

99

1

04/14/04

Scanned up to 370 nm


3 mW 99.2% Sc Rate

(500 nm) / 10 w/ 600 nm 3 mW

ethanol/acetone w/ other worked ok

Read and Understood By

Continued on Page



Signed

04/16/04

Date


Craig Chagge

Signed

4/21/04

Date

pattern issue number 4th character w/ Aragon
 on shore exposed by Santa Bay protect w/ methylene
 900m TMA + wood BOMA + 10m ABCN
 bubbled N_2 through for 15 min



1. 250/50 500 nm spacing 100 nm w/ 40 x 99% S

2. 250/25

Same

3. 50/10

Same

4. 500/250

washed Best

Same

500 nm spacing

Washed w/ ether

only see separation of features in #4

+ 10% BAPB in DMF for 15 min

Pattern came off glass. Should have

Coated 1st (washed to 50°C Blue N_2 Thresh)

Pattern DMF 3x

14 day more 10% BAPB + 500 nm DMF

5:50 am - 5:40 pm

Coated 30 DMF

Filled w/ Debye

Observed 2 systems on plate #4

□ □ □

DMF → □ □ □ DMF 10% SR 10 lines

□ □ □

10 day DMF + 10% BAPB + 500 nm DMF 18 min
 Coated 20x by DMF to get dye off (Continued on Page)

Read and Understood By

Signed

Date

Signed

Date

PROJECT _____

notebook No. W/ 100 MW


93

Continued From Page

370am

box no further

The two positions reported
are more fluorescent than
other positions though swimming
others are more fluorescent
(thick scutter) & none are
non fluorescent (in complete
protection).



2nd effect @ pattern by

or Notebook No. _____

Continued From Page _____

This approach:

- (1) repeat but cap w/ AA
- (2) make more hydrophobic polymer (increase from PEG w/ PEG)

Try both

TMX / PEG / GMA 45/45/10
 950m / 450mL / 45mL 11mg AIBN

~~DESS~~ 15m (BVB66L m₂)

used m₂
 10 Fm 5m

- (1) partial 500/25 w/ 1mm spaces @ 100mW 50% SR
- (2) wash excess w/ ether
- (3) blow N₂ through! Heat to 50°C for 15m
- (4) mab2
- (5) 10% BAPS w/ DMF 30m
- (6) N₂ 30m
- (7) AA 15m
- (8) repeat some pattern factors w/ 3m
- (9) Purulate 15m
- (10) Rise many times w/ DMF

1,4-Bis(3-oxo-2-propenyl)butane

Continued on Page _____

Read and Understood By


Signed

Date

Signed

Date

① The mes is getting 60 counts ~ 6 hrs what it was before
retention of the is from previous pre-hole or that the
pel is fluorescent.




①	500/25	w/ 100mm spray	792 SK
②	"	"	662 SK
③	"	"	332 SK
④	"	"	12 SK

#3 is more fluorescent (2x) than ①+②. Why? because of an
draining stream

② - wrong 5x fluorescent back: for the w/ E.P.H.
w/ E.T.A.Z.

bleed by w/ N₂
solid cluster

③ we moved to 5502 for 15 min while blowing
w/ through over



3000 w/
table @ 2000
0 370 m


Page

x/04

Date

PROJECT _____

① The mes of getting 60 counts ~ to this what it was before
notion of this is from previous Archival or that the
plate is fluorescent.



①	5500/25	w/ 1000nm spray	992 SR
②	"	"	662 SR
③	"	"	332 SR
④	"	"	172 SR

#3 is more fluorescent (2x) than ① & ② which are
becoming slower

② works 5x faster! for the w/ p.p.H
w/ CHAZ

bleed dry w/ N_2

Seal chamber in vacuum oven

③ worked to 55% for 15 min while blowing
 N_2 through oven

Revised 5x spray on - out

④ Revised ⑤ by MOC + 10ul AP24 + 500ul surf 30min on
table @ 200 rpm

WX objective 992 SR 1000nm input power @ 370 mW
features are ~ 9 counts

⑤ CAP: 146³⁰ 5ul AP24 + 20ul AA + 1ml surf 15 min
on orbital table @ 200 rpm

Continued on Page _____

Read and Understood By

Signed

Date

Signed

Date

⑦

Defect picked 4

50%


25%



10 lines for identification

LT 352

w/ 10 x @ 2 mm. net



12.10mm 4200/amp

992

1 1/2

Rough 300 in amp

⑧

10 day sample f. Sample after 1500 amp

15 min. on table @ 200 AMP

Consistent test for off of ketone. A larger

by night make better test w/ this system

Continued on Page

Read and Understood By

Signed

4/16/04

Date

Signed

4/21/04

Date

base to pattern by on photo polymer.


- I have been able to make small polymer features w/ the SS10 polymer, as small as 2 μ m. I have also made small features w/ my own methacrylate/AIBN formulations but I haven't quantified their sizes yet.

- I have imaged long 'hairs' of SS10 photo polymer w/ aspect ratios ~ 100 with a Confocal microscope. They appear to act as optical fiber like I see a large fluorescence when imaging at the glass polymer interface.

- I have successfully protected and deprotected my methacrylate resin functionalized w/ A diamine. The protection was w/ MOC. Deprotection w/ the same laser? Utilization using anhydride. This should provide significant improvements in my contrast ratio vs. working on functional glass. Further I hope the surface chemistry will be better vs. glass so I will be able to synthesize long peptides w/ high yield using the MOC protection/deprotection scheme.

Continued on Page

Read and Understood By



Signed

4/16/04

Date

Craig E. Mager

Signed

4/26/04

Date

purpose: to screen formulations against the following criteria

- ① rate of release of Absorbed dye (rate of washing)
- ② yield of MOC protected acetone

Methodology: prepare samples, split onto 96 well plates, absorbent to 50% for 15 min with blowing N₂ through over. React w/ Diamine, MOC, & ready to use on plate 94-95. React w/ solvent (DMF) & look for maximum fluorescence. Retest Best formulation for activity against Diamine.

Samples


All 12 ABN samples. make STD solution of 1% TMA/GMA

Trimethylolpropane (TMA) Glycidyl methacrylate (GMA)

#		
1	10	90
2	20	80
3	30	70
4	40	60
5	50	50
6	60	40
7	70	30
8	80	20
9	90	10
10	95	5
11	100	0

Continued on Page 99

Read and Understood By



Signed

4/20/04

Date



Signed

4/21/04

Date

PROJECT

#	Trimethylamine ethoxide (14/3 Eo/04)	Amr
	For-ethylate	
A	90	10
B	75	25
C	50	50
D	25	75
E	10	90
F	0	0

So marked each fr, & see
 soaked into each. Under placed in 70°C oven
 in the Blaney Throsh, in Glass Tray 11:09 AM

Added 10ul of 1% ~~GA~~ PB to each bucket 15 ml
all samples have swirled

A, B, & C Swelled the nest
223 maybe the best or the showiest.

All become white when swollen.

Wish to try more cross-linking. or it curls so
that I have not cond. it long enough or
don't purge w/ N_2 well enough.

was to do a series of the vast Aegle
and have a whole team of

Name		In Charge		Time		Amn	GMA	TMA
Time	Direct	Indirect	Tr. with	Tr. with	Tr. with			
High 3hr	90	90	90	90	90	60	10	
Low 15min	90	90	90	90	90	60	10	
Read and Understood By						10		
4/20/04						10		
Signed						10		
Date						10		

Continued on Page

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/015764

International filing date: 06 May 2005 (06.05.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US
Number: 60/608,774
Filing date: 10 September 2004 (10.09.2004)

Date of receipt at the International Bureau: 20 June 2005 (20.06.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1331929

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

June 09, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

APPLICATION NUMBER: 60/608,774

FILING DATE: *September 10, 2004*

RELATED PCT APPLICATION NUMBER: *PCT/US05/15764*



Certified by

Don W. Dudas

Under Secretary of Commerce
for Intellectual Property
and Director of the United States
Patent and Trademark Office

091004

22713
U.S. PTO

Express Mail Label No. EL988555789US

Please type a plus sign (+) inside this box → ☐Approved for use through 7/31/2003.
U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

17497 U.S. PTO
60/608774

091004


PROVISIONAL APPLICATION FOR PATENT COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

INVENTOR(S)					
Given Name (first and middle [if any])	Family Name or Surname	Residence (City and either State or Foreign Country)			
Trent Russell	Northen	Tempe, Arizona			
Additional inventors are being named on the ____ separately numbered sheets attached hereto					
TITLE OF THE INVENTION (500 characters max)					
LIGHT ACTIVATED MOVING POLYMER					
Direct all correspondence to: CORRESPONDENCE ADDRESS					
<input checked="" type="checkbox"/> Customer Number <input type="text" value="26707"/>					
OR Type Customer Number here					
<input type="checkbox"/> Firm or Individual Name					
Address					
Address					
City		State		ZIP	
Country		Telephone		Fax	
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification Number of Pages <input type="text" value="12"/>		<input type="checkbox"/> CD(s), Number <input type="text"/>		<input type="checkbox"/> Other (specify) <input type="text" value="Cover sheet; Postcard"/>	
<input type="checkbox"/> Drawing(s) Number of Sheets <input type="text"/>					
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76					
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT					
<input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.					FILING FEE AMOUNT (\$)
<input type="checkbox"/> A check or money order is enclosed to cover the filing fees					<input type="text" value="\$80"/>
<input checked="" type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: <input type="text" value="17-0055"/>					
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.					
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.					
<input checked="" type="checkbox"/> No.					
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____					

Respectfully submitted,

SIGNATURE



TYPED or PRINTED NAME

Robert D. Atkins

TELEPHONE (602) 229-5311

Date 09/10/04

REGISTRATION NO.
(if appropriate)
Docket Number:

34,288

112624.00147

USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is used by the public to file (and by the PTO to process) a provisional application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the complete provisional application to the PTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Patent Application, Commissioner for Patents, Alexandria, VA 22313-1450.

1870507

EXPRESS MAIL LABEL NO.: EL988555789US

I hereby certify that this correspondence listed below is being deposited with the United States Postal Service on the date set forth below as Express Mail in an envelope addressed to: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

Date of Signature
and Deposit: September 10, 2004

By:

(Signature of person depositing mail)

MARITZA O'NEILL

CERTIFICATE OF MAILING PURSUANT TO 37 C.F.R. 1.10

Applicant: Trent Russell Northen

Date of Filing: September 10, 2004

Title: *LIGHT ACTIVATED MOVING
POLYMER*

Art Unit: Unassigned

Examiner: Unassigned

Attorney Docket No.: 112624.00147

Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

Type of Filing:

- 1) Provisional Application For Patent Cover Sheet
- 2) Specification (12 pages, plus cover sheet)
- 3) Return postcard

PATENT

PROVISIONAL APPLICATION

Of

TRENT RUSSELL NORTHEN

For

UNITED STATES LETTERS PATENT

on

LIGHT ACTIVATED MOVING POLYMER

Attorneys:

QUARLES & BRADY STREICH LANG L.L.P.

ONE RENAISSANCE SQUARE

TWO NORTH CENTRAL AVENUE

PHOENIX, AZ 85004-2391

Express Mail Label No.: EL988555789US
Attorney Docket No.: 112624.00147

Light Directed Movement of Polymer Microstructures

Trent R. Northen, Neal W. Woodbury,

Department of Chemistry and Biochemistry, Arizona State University, Biodesign Institute at Arizona State University, Tempe, AZ, 85283.

Abstract

Light induced surface chemistry changes have been used to move swollen polymer microstructures. Swellable trimethylolpropane trimethacrylate (TRIM) crosslinked poly(2-hydroxyethyl methacrylate) conical microstructures were constructed by azo-bis-isobutyronitrile (AIBN) photopolymerization using a 20x 0.5NA microscope objective and 365nm laser excitation. Structures were aminated with glycine and protected with the photolabile group 4-nitroveratryloxycarbonyl (NVOC). Differential swelling with and without NVOC of 10% was observed in N,N'-dimethylformamide (DMF). Removal of NVOC with 365nm laser excitation induced polymer shrinkage in excess of 4%, resulting in maximum polymer velocity of 1mm/s, and displaced solvent velocities in excess of 0.01mm/s.

Introduction

A number of synthetic polymers have recently been developed that respond to changes in surface energy resulting from external stimuli including: mechanical deformation, heating, solvent contact, and exposure to light[1]. Exciting applications of such materials include: implants based on shape-memory materials, gels respond (e.g. swell) in response to changes in pH or specific molecules may be used for feedback control for drug delivery, and microfabricated vascular networks[2]. Given the fact that three-dimensional polymer structures can now be constructed on the submicron scale using nonlinear laser patterning [3, 4] [5, 6], it should also be possible to develop micro or nanomechanical devices based on polymer movement.

A particularly versatile stimulus that could be used for directing polymer movement at dimensions down to the submicron level is light. Photolabile protective groups offer the ability to selectively break bonds using light and therefore substantially change the surface characteristics of the polymer in a light-directed fashion. 4-nitroveratryloxycarbonyl (NVOC) is a common photolabile group and is known to cleave using a Norrish-type II reaction [7]. It has found wide use in protecting amines [8] and has applications including: photogeneration of organic bases [9], microarrays [10], novel proteins [11], and variations of NVOC as linkers in peptide synthesis [12]. In these cases the addition and removal of NVOC modulates the reactivity of an amine. This work describes the use of NVOC to instantaneously modulate the surface properties of a porous polymer using light.

Porous polymers are common in solid phase synthesis [13], drug delivery [14, 15], tissue engineering [16], and separations [17] [18]. A range of polymers are now used for solid phase synthesis including polyacrylate resins [19]. These porous polymer structures are the result of phase separation during free radical crosslinking copolymerization swell in

Light Directed Movement of Polymer Microstructures

compatible solvents [20]. Typically solvents (porogens) can be used to control the pore size [21]. The surface of the polymer is often modified to improve functionality [22] or as a result of solid phase synthesis.

Here we report polymer microstructures that shrink and move when illuminated with light. Cleavage of NVOC immediately exposes the primary amine, resulting in large changes in the surface chemistry and swelling of the polymer. This allows light-directed spatial control of polymer movement.

Experimental Details

Materials: Glass coverslips for an FCSII chamber (see below) were purchased from Biopetechs (Butler, PA). 2-hydroxyethyl methacrylate (HEMA), trimethylolpropane trimethacrylate (TRIM), azo-bis-isobutyronitrile (AIBN), piperidine, diisopropylethylamine were purchased from Sigma-Aldrich Chemical Co. (Milwaukee, WI). 4-nitroveratryloxycarbonyl chloride and 3-(trimethoxysilyl)propyl methacrylate were from Fluka GmbH (Buchs, Switzerland). Dimethylformamide (DMF) was from Applied Biosystems Inc. (Foster City, CA). Methanol, hydrogen peroxide (30%), sulfuric acid, hydrochloric acid were purchased from Mallinckrodt Inc. (Paris, KY). Isopropanol and ethanol were from ACROS Organics (Geel, Belgium). Acetonitrile was from Alfa Aesar (Ward Hill, MA). 2-(1H-benzotriazole-1-yl)-1,1,3,3-tetramethyluronium hexafluorophosphate (HBTU) and Fmoc-glycine (Fmoc-Gly) were from Advanced ChemTech Inc. (Louisville, KY). Finally, water was purified using a NANOPure ultrapure filtration system Barnstead. (Dubuque, IA).

Equipment: All reactions were performed inside an FCSII flow chamber Biopetechs Inc. (Butler, PA). Patterning and deprotection were done using light from a mode-locked Tsunami Ti:sapphire laser pumped by a 5 W Millennia Vs diode-pumped cw laser, Spectra-Physics Inc. (Mountain View, CA), through a 20x 0.5NA objective attached to an Eclipse TE2000-U microscope, Nikon Inc. (Japan) equipped with a ProScan microscope stage, Prior Scientific Inc. (Rockland, MA). The laser beam was modulated using a Model 350-80 electro-optic light modulator with model 302 power supply, Conoptics Inc. (Danbury, CT) controlled by software developed in-house. Laser power was measured using a Model 1815-C power meter, Newport Co. (Irvine, CA). Images taken using Cascade Photometrics CCD, Roper Scientific Inc. (Tucson, AZ) through 10x 0.3NA objective lens, Nikon Inc. (Japan) using MetaVue 6.0 software, Universal Imaging Corporation Limited (Marlow, UK) for acquisition and analysis. Scanning electron microscopy (SEM) was performed using a XL30ESEM environmental SEM, FEI Co. (Hillsboro, OR) on a sample coated with 3.5nm palladium/gold.

Surface Functionalization. Glass cover slides for a FCSII flow chamber were cleaned using a modification of methods reported by McGall [23]. Briefly: slides were soaked 15 min at RT with 60/40 (v/v) sulfuric acid/hydrogen peroxide (use extreme caution when using this solution), placed in 10% sodium hydroxide (w/v) at 70°C for 3 min and placed in 1% HCl at RT for 1 min. Between each step the slides were soaked in nanopure water for 3 minutes. A solution of 1% 3-(trimethoxysilyl)propyl methacrylate in 95% ethanol 5 % water was mixed for 10 minutes, and the slides were reacted at RT for 15 minutes with gentle agitation. Slides were then soaked in isopropyl alcohol for 3 min, nanopure water for 1 min, and then placed in a 100°C oven for 5 minutes after which the oven was turned off and nitrogen was blown through for 1 hr. The slides were stored under nitrogen until they were used.

Light Directed Movement of Polymer Microstructures

Fabrication of Polymer Structures. A total of 6mg of AIBN was dissolved in 95 μ L HEMA and 579 μ L TRIM. This was placed in an optical chamber, and irradiated with 4 mW (all powers reported are measured entering microscope) of 365nm (8nm full-width-at-half-maximum) light for 1.6s per feature through a 20x objective focused 400 μ m above the surface of the cover slip. Excess monomer was drained and sample washed with methanol and DMF. The features were spaced 600 μ m apart..

Amination of Microstructures. Fmoc-Gly was coupled to the photopolymer hydroxyl group using 18.6mg Fmoc-Gly, 22.5mg HBTU, 11.5 μ L DIPEA, and 600ul DMF.

Reaction was mixed at 50°C for 30 min. The structures were then rinsed with DMF and the Fmoc removed with 20% piperidine in DMF for 10min. The yield of the reaction was determined using the absorbance at 301nm for the Fmoc-piperidine adduct. Typical polymer substitution levels were 0.1 nanomoles/feature.

Coupling NVOC and 6-nitrophenyl chloroformate (NPC) to aminated microstructures. A solution of 19mg NVOC or 14mg NPC, 40 μ L DIPEA, and 600 μ L DMF was reacted with polymer microstructures by mixing for 30min at 50°C.

Laser cleavage of NVOC: The same laser beam used for making the microstructures was used for cleavage of the NVOC. The beam was attenuated as needed.

Swelling and Tip Velocity Measurements: Images taken of the microstructures at the glass polymer interface in various solvents were manually fitted with ellipses of known pixel area. Tip velocity was calculated from the distance moved in sequential images over a known amount of time.

Results and Discussion

Porous polymer microstructures. Polymer structures were obtained via the photopolymerization of HEMA and TRIM with AIBN. Oxygen quenching[6] and light
Light Directed Movement of Polymer Microstructures

attenuation from AIBN absorption were used to limit the polymer structure dimensions to the volume of excitation between the surface and the focus of the laser.




Figure 1. SEM images of polymer microstructures: (A) a portion of the array, (B) one microstructure, (C) top of microstructure, (D) macropore at tip of microstructure, (E) macropores below tip, (F) junction of smooth central region with lowest rough/microporous region.

An array of polymer features was generated by laser-directed photopolymerization (Figure 1A). Partial polymerization results in soft porous structures that were measured on the optical microscope to be 400 μm tall, having an elliptical base with radii of 75 and 200 μm . The structure of each feature has a heterogeneous morphology due to spatial differences in light intensity in the focused laser beam. Structures appear composed of four regions (Figure 1B): two macroporous regions close to the beam focus (Figure 1C-E), an apparently nonporous central region (smooth area in Figure 1F), and a rough potentially porous region nearest the glass surface (rough area in Figure 1F). The macropores at the top are on order of 1 μm . The number of reactive sites (using FMOC as a probe) was estimated to be 0.1 nmole per feature, 5 orders of magnitude more than would be expected for a nonporous material.

Solvent Swelling: The swelling of the polymer with and without the NVOC protective group was measured in various solvents. It was found that the greatest swelling, and largest difference in swelling between the protected and unprotected resin, was in the polar aprotic solvents DMF and acetonitrile (Figure 2). The swelling of the polymer was found, as expected, to be related to

Light Directed Movement of Polymer Microstructures

BEST AVAILABLE COPY

the Hildebrand solubility parameter [24] where maximum swelling of a slightly crosslinked polymer occurs in solvents of similar solubility parameter [24].




Figure 2. (top) Images of NVOC-Gly polymer features laying on the glass surface, illustrating the difference in swelling when exposed good solvents (acetonitrile) and poor solvents (methanol). The images were recorded using a 10x objective lens focused at polymer glass
Light Directed Movement of Polymer Microstructures

interface. (bottom) swelling of polymer, Glycine-polymer, NVOC-Gly polymer in various solvents as a function of the Hildebrand parameter.

It is apparent from Figure 2 that the swelling of the resin changes dramatically with the solvent. In DMF the area of the NVOC-Gly polymer at the glass polymer interface increased 10% versus the Gly polymer. This has been seen before; it was found by Merrifield that over the course of solid phase peptide synthesis, resin swelling increased more than five fold [25]. This behavior was attributed to the net decrease in free energy upon swelling due to solvation of peptide chains bound to the polymer matrix. Presumably a similar solvation mechanism accounts for the differential swelling of the resin with and without NVOC as shown in figure 2.

Polymer Movement. Upon laser excitation of the NVOC-Gly polymer structures in DMF or acetonitrile, the NVOC is photocleaved, resulting in shrinking of the illuminated portion of the polymer, causing the polymer to bend. Figure 3 shows a series of images collected as a polymer structure moves towards the laser beam.





Figure 3. Progression of images of NVOC-Gly polymer structure with asymmetric illumination at lower edge of polymer structure. Perspective is looking down on the polymer as it bends toward the lower edge of the images. The upper left image of the series represents the point at which illumination was initiated, with time continuing from left to right and then down the array. Images were taken through a 10x objective lens at 5ms per frame and the solvent used was acetonitrile.

Light Directed Movement of Polymer Microstructures

BEST AVAILABLE COPY

Symmetric illumination of an NVOC-Gly polymer feature results in the rapid release of solvent from the microstructure. Small particulates in solution rapidly move radially away from the polymer structure during shrinkage with a maximum velocity of 0.01 mm/s. This provides a lower limit for the maximum fluid velocity. It should also be noted that the NVOC group cleaves as nitrosobenzyl aldehyde [Patchornik, 1970 #1] which is released upon illumination. A modified photocleavable group of this nature could be used as a method for local delivery of reagents or drugs.

The base area of the polymer shrinks by ~4% after a 20sec 400uW illumination period (Figure 4 top and bottom respectively). Given that the polymer is covalently attached to the glass it is expected that the actual bulk shrinkage is greater than is observed at or near the glass surface.



Light Directed Movement of Polymer Microstructures

Figure 4. (top) Movement of particles away from polymer upon illumination with light.

(bottom) Difference image showing shrinkage of polymer structure upon laser excitation, before (black) and after (white) 20sec of 400 μ W

The polymer movement is very rapid especially at the tip of the polymer structure. Velocities on order of 1mm/sec were recorded (Figure 5) in acetonitrile. This is several times faster than that observed in DMF (Figure 5) even though the differential swelling is larger in DMF (Figure 1). This is explained by the three fold higher viscosity of DMF resulting in greater resistance to flow, slowing the movement of the microstructure. It is also clear that there was little if any movement of structures in solvents that did not swell the resin (Figure 5).




Figure 5. Tip velocity as a function of excitation power in various solvents.

Chemical mechanism of polymer volume change upon illumination: The design hypothesis underlying the development of this system was that photocleavage induced changes in surface chemistry result in the volume changes and associated polymer movements. The other possible mechanisms for mechanical movement upon illumination and photocleavage include photochemical curing processes, electrostatic forces, hydrogen bonding, and optical trapping, and measurements were performed to investigate each of these.

A photochemical curing process is unlikely for several reasons. The polymer movement continues for ~1sec after a 100ms 1mW exposure in DMF. Given the presence of oxygen and low viscosities of the solvents it is unlikely that any free radical or photochemical process would continue for this long in the dark. It was also found (Figure 6) that these dramatic movements

Light Directed Movement of Polymer Microstructures

only occurred with NVOC and not in the presence of the polymer itself or with another chromophore (NPC). It was possible to reattach NVOC and regain partial polymer movement.




Figure 6. Tip velocity as a function of input power for polymer structures in DMF.

To test for the role of hydrogen bonding of the hydroxyl or amine groups, the polymer was soaked for approximately 1 hour in acetic anhydride and then illuminated. While this reduced the speed of tip movement (50 $\mu\text{m/s}$ with 400 μW excitation), the movement was still substantial, showing that the movement is not a result of hydrogen bond formation upon NVOC cleavage since all hydroxyls and amines should have been rapidly acylated in acetic anhydride.

The fact that the polymer curves towards the laser beam makes an electrostatic mechanism (formation of protonated amine groups or charged intermediates in the photocleavage reaction) unlikely since one would expect repulsive electrostatic forces to push the polymer away from the illuminated region (Figure 3). Additionally, the movement was not observed in nonpolar solvents, where the electrostatic effect should be greatest and was observed even in a solvent system with a 60 mM ionic strength that should strongly shield the charge-charge interactions (data not shown). Finally, the movement of the polymer is not dependent on the position of the focus, making optical trapping a very unlikely mechanism. Polymer movement is observed even when the focus is at the polymer glass interface, 400 μm from the polymer tip.

Conclusion

We have described porous polymer structures that make dramatic movements with the rapid release of solvent when illuminated. We hope to attach spiropyrans to these polymer structures in attempt to make a reversible system. This may find use in systems where it is desirable to control the movement of a polymer structure or for releasing material into solution. Such a system would be very exciting for converting light energy into mechanical movement.

Light Directed Movement of Polymer Microstructures

Acknowledgements

Sudhir Gudala for developing the scanning software. Saritha Keshammolu for surface preparation.

Center for Solid State Science at Arizona State University for SEM

References

1. Russell, T.P., *Surface-responsive materials*. Science, 2002. 297(5583): p. 964-967.
2. Langer, R. and D.A. Tirrell, *Designing materials for biology and medicine*. Nature, 2004. 428(6982): p. 487-492.
3. Kawata, S., et al., *Finer features for functional microdevices - Micromachines can be created with higher resolution using two-photon absorption*. Nature, 2001. 412(6848): p. 697-698.
4. Sun, H.B., et al., *Three-dimensional nanonetwork assembled in a photopolymerized rod array*. Advanced Materials, 2003. 15(23): p. 2011-2014.
5. Belfield, K.D., et al., *Multiphoton-absorbing organic materials for microfabrication, emerging optical applications and non-destructive three-dimensional imaging*. Journal of Physical Organic Chemistry, 2000. 13(12): p. 837-849.
6. Maruo, S. and K. Ikuta, *Three-dimensional microfabrication by use of single-photon-absorbed polymerization*. Applied Physics Letters, 2000. 76(19): p. 2656-2658.
7. Bochet, C.G., *Photolabile protecting groups and linkers*. Journal of the Chemical Society-Perkin Transactions 1, 2002(2): p. 125-142.
8. Patchornik, A., B. Amit, and R.B. Woodward, *Photosensitive protecting groups*. Journal of the American Chemical Society, 1970. 92(21): p. 6333-5.
9. Cameron, J.F. and J.M.J. Frechet, *Photogeneration of Organic-Bases from Ortho-Nitrobenzyl-Derived Carbamates*. Journal of the American Chemical Society, 1991. 113(11): p. 4303-4313.
10. Fodor, S.P.A., et al., *Light-Directed, Spatially Addressable Parallel Chemical-Synthesis*. Science, 1991. 251(4995): p. 767-773.
11. Mendel, D., J.A. Ellman, and P.G. Schultz, *Construction of a Light-Activated Protein by Unnatural Amino-Acid Mutagenesis*. Journal of the American Chemical Society, 1991. 113(7): p. 2758-2760.
12. Holmes, C.P. and D.G. Jones, *Reagents for Combinatorial Organic-Synthesis - Development of a New O-Nitrobenzyl Photolabile Linker for Solid-Phase Synthesis*. Journal of Organic Chemistry, 1995. 60(8): p. 2318-2319.
13. Merrifield, R.B., *Solid Phase Peptide Synthesis .I. Synthesis of a Tetrapeptide*. Journal of the American Chemical Society, 1963. 85(14): p. 2149-&.
14. Kost, J. and R. Langer, *Responsive Polymer Systems for Controlled Delivery of Therapeutics*. Trends in Biotechnology, 1992. 10(4): p. 127-131.
15. Langer, R., *Drug delivery and targeting*. Nature, 1998. 392(6679): p. 5-10.
16. Liu, V.A. and S.N. Bhatia, *Three-dimensional photopatterning of hydrogels containing living cells*. Biomedical Microdevices, 2002. 4(4): p. 257-266.
17. Svec, F., et al., *Design of the monolithic polymers used in capillary electrochromatography columns*. Journal of Chromatography A, 2000. 887(1-2): p. 3-29.
18. Righetti, P.G., *Macroporous Gels - Facts and Misfacts*. Journal of Chromatography A, 1995. 698(1-2): p. 3-17.
19. Kempe, M. and G. Barany, *CLEAR: A novel family of highly cross-linked polymeric supports for solid-phase peptide synthesis*. Journal of the American Chemical Society, 1996. 118(30): p. 7083-7093.

Light Directed Movement of Polymer Microstructures

20. Okay, O., *Macroporous copolymer networks*. Progress in Polymer Science, 2000. 25(6): p. 711-779.
21. Yu, C., et al., *Preparation of monolithic polymers with controlled porous properties for microfluidic chip applications using photoinitiated free-radical polymerization*. Journal of Polymer Science Part a-Polymer Chemistry, 2002. 40(6): p. 755-769.
22. Rohr, T., et al., *Photografting and the control of surface chemistry in three-dimensional porous polymer monoliths*. Macromolecules, 2003. 36(5): p. 1677-1684.
23. McGall, G.H., et al., *The efficiency of light-directed synthesis of DNA arrays on glass substrates*. Journal of the American Chemical Society, 1997. 119(22): p. 5081-5090.
24. Barton, A.F.M., *Solubility Parameters*. Chemical Reviews, 1975. 75(6): p. 731-753.
25. Sarin, V.K., S.B.H. Kent, and R.B. Merrifield, *Properties of Swollen Polymer Networks - Solvation and Swelling of Peptide-Containing Resins in Solid-Phase Peptide-Synthesis*. Journal of the American Chemical Society, 1980. 102(17): p. 5463-5470.

This invention can be used as _ This technology would allow the conversion of light energy to mechanical energy either through the movement of a fluid or the polymer itself. There would seem to be a wide range of possibilities for such a material. Ranging from large solar collectors, light powered nanobots (photopolymer structures 100's of nanometers can be readily made), artificial muscle, drug delivery systems, microfluidic pumps and valves, etc.

1. This invention provides the following advantages: It results in dramatic changes in the physical dimensions of the polymer.
2. It releases (or could absorb) solvent
3. It is a general system that could be used with any porous polymer formulation or potentially on the surface of very thin polymer structures.
4. The surface area of a porous polymer is many orders of magnitude higher than a non-porous polymer, this would be the preferred mode for drug delivery

There are molecules (Azobenzene, spiropyrans, etc) that act as molecular switches, one color of light puts them in one form, another moves them back to the initial form. By attaching one of these molecules to the surface that has a large polarity change upon switching forms, it should be possible to make a polymer that expands with one color of light and contracts with another color of light. This technology would allow the conversion of light energy to mechanical energy either through the movement of a fluid or the polymer itself. There would seem to be a wide range of possibilities for such a material. Ranging macroscale solar collectors, light powered nanobots (photopolymer structures 100's of nanometers can be readily made), artificial muscle, drug delivery systems, microfluidic pumps and valves, etc.

There is a NVOC derivative that is used as a linker in peptide synthesis. It could be used to release a material of interest. One end of the linker would be linked to the polymer, the other to the material (drug) to be released. By adjusting the surface energy of the polymer it would be possible to design a system that would rapidly release the material with light. Merrifield has shown that polymer resin swells 5x with a large peptide attached—this system with a photocleavable linker would allow the delivery of peptide (and other) drugs accompanied with a rapid movement of solvent.

It may be possible to use a conducting polymer and switch the polymer states by oxidizing and reducing groups on the surface electrically. This would allow this technology to be used in places not accessible to light (inside the body) or in electrical devices.

There are other photoactivated groups and polymers that could be used.

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/015764

International filing date: 06 May 2005 (06.05.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US
Number: 60/623,181
Filing date: 29 October 2004 (29.10.2004)

Date of receipt at the International Bureau: 20 June 2005 (20.06.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1331930

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

June 09, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

APPLICATION NUMBER: 60/623,181

FILING DATE: *October 29, 2004*

RELATED PCT APPLICATION NUMBER: *PCT/US05/15764*



Certified by

Under Secretary of Commerce
for Intellectual Property
and Director of the United States
Patent and Trademark Office

102904

16138 U.S. PTO

Approved for use through 07/31/2006. OMB 0651-0032
U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

PROVISIONAL APPLICATION FOR PATENT COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. EL988555877US


(09/04)

19249 U.S. PTO
60/623181

102904

INVENTOR(S)		
Given Name (first and middle [if any])	Family Name or Surname	Residence (City and either State or Foreign Country)
Trent Russell Neal Walter	Northen Woodbury	Tempe, Arizona Tempe, Arizona
Additional inventors are being named on the _____ separately numbered sheets attached hereto		
TITLE OF THE INVENTION (500 characters max):		
PEPTIDE CHARACTERIZED FOR PATTERNED PHOTOPOLYMER FORMED USING LIGHT DIRECTED SYNTHESIS		
Direct all correspondence to: CORRESPONDENCE ADDRESS		
<input checked="" type="checkbox"/> The address corresponding to Customer Number: 26707		
OR		
<input type="checkbox"/> Firm or Individual Name		
Address		
City	State	Zip
Country	Telephone	Fax
ENCLOSED APPLICATION PARTS (check all that apply)		
<input checked="" type="checkbox"/> Specification Number of Pages 21		
<input type="checkbox"/> Drawing(s) Number of Sheets _____		
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76		
<input type="checkbox"/> CD(s), Number of CDs _____		
<input checked="" type="checkbox"/> Other (specify) Cover Sheet; postcard		
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT		
<input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.		
<input type="checkbox"/> A check or money order is enclosed to cover the filing fees.		
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.		
<input checked="" type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: 17-0055		
A duplicative copy of this form is enclosed for fee processing.		
<input checked="" type="checkbox"/> invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.		
No.		
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____		

SIGNATURE



TYPED or PRINTED NAME Robert D. Atkins

TELEPHONE 602-229-5311

Date October 29, 2004

REGISTRATION NO. 34,288

(if appropriate)

Docket Number: 112624.00138

USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.11 and 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

1882295

EXPRESS MAIL CERTIFICATE: EL988555877US

I hereby certify that this correspondence listed below is being deposited with the United States Postal Service on the date set forth below as Express Mail in an envelope addressed to: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

Date of Signature
and Deposit: October 29, 2004

By: 

(Signature of person depositing mail)

MARITZA O'NEILL

CERTIFICATE OF MAILING PURSUANT TO 37 C.F.R. 1.10

Applicant: Trent R. Northen et al.

Date of Filing: October 29, 2004

Title: PEPTIDE CHARACTERIZED FOR
PATTERNED PHOTOPOLYMER
FORMED USING LIGHT DIRECTED
SYNTHESIS

Art Unit: Unassigned

Examiner: Unassigned

Attorney Docket No.: 112624.00138

Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

Type of Filing:

- 1) Provisional Application For Patent Cover Sheet
- 2) Specification (21 pages, plus cover sheet)
- 3) Return postcard

PATENT

PROVISIONAL APPLICATION

of

TRENT RUSSELL NORTHEN

NEAL WALTER WOODBURY

For

UNITED STATES LETTERS PATENT

on

**PEPTIDE CHARACTERIZED FOR PATTERNED PHOTOPOLYMER FORMED
USING LIGHT DIRECTED SYNTHESIS**

Attorneys:

QUARLES & BRADY STREICH LANG L.L.P.

ONE RENAISSANCE SQUARE

TWO NORTH CENTRAL AVENUE

PHOENIX, AZ 85004-2391

Express Mail Label No.: EL988555877US
Attorney Docket No.: 112624.00138

BEST AVAILABLE COPY

Prep: to pattern Fmoc-GGFLC(OMe) on the polymer
that represents w/ Muc & pattern
P.S.U.
Hydroxide Antibody & detect via fluorescence.

From page 29 9.5 mL H₂O + 5.72 mL DMF + 4 mg A.B.N.
Used & Sonicated until dissolved
added to Chamber w/ Rough methacrylate slide
Patterned at 200 730/2 8.0 mFWHM 8.0 W 2 sec 525 nm
2x 27x13 arrays, looks great!

Acidolysis Fmoc-GGFLC(OMe) prep'd by A.B.N. still
dissolved in TFA/THF/succinic anhydride 1:2:1 = 80 mg
extracted w/ MeOH (aq. KOH w/ 60% D.P.E.A. in
MeOH unnecessary) Dried down from layer (yellow)
Bong - 60 mg (7.42 mg) = 35.5 mg peptide
Fmoc - G - G - F - L (OMe)

$$\text{Fmoc } 223.3 + 2(57.05) + 147.18 + 113.16 + 17 = 614.74 \frac{\text{g}}{\text{mole}}$$

$$35.5 \text{ mg} / 614.74 \frac{\text{g}}{\text{mole}} = 0.058 \text{ mmole} \times \left(\frac{2.35}{2.5} \right) = 5.5 \times 10^{-2} \text{ mmole D.P.E.A.}$$

$$\times 28 = 0.29 \text{ mmole D.P.E.A.}$$

$$(5.5 \times 10^{-2} \text{ mmole}) (379.25 \frac{\text{mg}}{\text{mole}}) = 20.9 \text{ mg HSTU}$$

$$0.29 \text{ mmole} (127.25 \frac{\text{mg}}{\text{mole}}) (1.4 \frac{\text{mL}}{\text{mg}}) = 52.5 \mu\text{L}$$

9.11.04 9.15.04

30 mM	614.74 mg/mole	= 9.22 mg	9.5 mg p-peptide
-------	----------------	-----------	------------------

$$(30 \text{ mM}) (379.25 \frac{\text{mg}}{\text{mole}}) = 5.68 \text{ mg HSTU}$$

$$(60 \text{ mM}) (127.25 \frac{\text{mg}}{\text{mole}}) (6.742 \text{ mg}) = 5.2 \mu\text{L D.P.E.A.}$$

Use 150 mM MeOH Read and Understood By
9.15.04 10/29/04
Signed _____ Date _____ Signed _____ Date _____

800

720

800

Mass (mg)

Hot

Continued on Page 72

Sample B, Botaine-GGFL-omide
D.L. 10/29/04


Followed procedure from T-8-71

Patterning looks great

- (1) Resist 30 mils on PMF
- (2) Resist/HSTU/PMF 11:00 am - 12:09 pm 50°C 150 rpm
- (3) 1hr @ 50°C PMF 3x PMF 10mm PMF $OD_{301} = 0.053$
- (4) 10mm PMF $OD_{301} = 4 \times 1.24 = 4.96$
- (5) Resist 4x PMF
- (6) Fmoc Y66 EL Coat ~ 15µ / 13µ / 1.5µ HSTU / 50ml PMF
- (7) 2pm 9.29.04 50°C 150 rpm
- (8) 4pm Resist 2x @ 50°C 2x (7x PMF)
- (9) 10mm Resist $OD_{301} = 0.078$
- (10) 10mm PMF $OD_{301} = 1.7$
- (11) 3x PMF
- (12) 14ms NDC 40ml Resist 600ml PMF 6pm 150rpm 50°C
- (13) Resist 3x PMF 7pm left in PMF overnight
- (14) 830 pm 1/30 Resist 3x PMF 3x ACN
- (15) Added 600ml ACN + 20ml D₂O
- (16) Wash as before 730/2 30m Resist
- (17) adjust to 2ml water w/ 200 20 sec exposure
- (18) See lower setting
- (19) It was observed that the Resist no longer seemed
- (20) Wash in Acetonitrile, w/ 200 let more in buffer
- (21) to replace w/ ACN 70823
- (22) Washed 3x w/ PMF ACN PMF ACN = 150471
- (23) Wash 2x reference
- (24) Resist 3x ~~PMF~~ ACN 3x PMF, 3x ACN
- (25) Dried w/ N₂

Continued on Page _____

Read and Understood By



Signed

10/27/04

Date




Signed

10/29/04


Date




Voyager Spec #1=>MC[BP = 964.4, 35809]



ISO: C48H63N5O14P



TMPP + GGFL papillo
D:\VOSTNref_0007.dat
Acquired: 12:40:00, October 05, 2004




GRGDSP

PROJECT

Notebook No. _____

85

Continued From Page _____



(1.4 mg) (100 mL) (76.8 mg) = 0.0768 mg TMP-AC - OBT BR -
 from P. Brown Ref. Analytical Biochem
 68 305-317 1999

5x = 0.384 mg

~ 0.5 mg TMP-AC - OBT BR - + 20 mL DMSO + 480 mL DMF
 10 mL stored @ 35°C 150 RPM

(2) 11 to 20 mL rinsed 30 mL left @ 150 RPM 35°C - 1 pm
 rinsed 3x more! Rinsed 1/2 of each array
 placed in 2 separate w/ 100% 2 Tri: 100% 2 S. line / TFA
 Sonicated for 15 min light exposed send to EIR
 apart more than in exposed.

Concentrated! Came to hand

Mass spec clearly shows the TMP-66FL mass $C_{43}H_{63}N_5O_{14}P$ 946.401
 there are the correct 110.02 peaks! The post source decay
 shows TMP-66FL, TMP-66F, TMP-66, TMP-6
 - this was successful. There was a small amount of TMP-66FL
 in the Control, however. The major peak was the 653.24 1.04
 corresponding to mass-66FL. This demonstrates the
 successful light directed synthesis of a polymer on a photo-
 patterned polymer support & the mass spec characterized
 of the product from the polymer material. It may be
 possible to do the mass spec in situ by spraying matrix/TFA
 directly on the plate (control) - leaving into mass spec

Read and Understood By

Signed

10.27.04

Date

Signed

10/29/04

Date

9E4

30821

PROJECT _____


Notebook No. _____

Continued From Page _____

72⁷³

9/20/2004 1:25:05 PM

Page 1 of 24



10min rinse2
 10min piperidine
 10min rinse after
 rinse after coupl
 10min piperidine
 10 min rinse afte
 10 min deprotect
 10 min rinse afte
 10 min rinse afte
 10 min piperidin

rt

72m
 600
 600
 550
 450
 420
 360
 300
 240
 180
 120
 60
 0

Signed _____

Date _____

Signed _____

Date 10/29/04

Process: to pattern Fmoc-66FL CooH on the polymer
that represents w/ MSC & pattern
P.S.U.
Hydrazine Antibody & detect the fluorescence.

From page 29 95% Hexam + 5% 2nd trim + 4mg ABN
Used & Sonicated until dissolved
added to chamber w/ growth medium w/ slide
patterned on page 29 200 730/2 800FWHM 80W 2500nm
2x 27x13 arrays, looks great!

to use Fmoc-66FL CooH prepared by DAN. Still
dissolved in TFA/THF:supernatant solution. 1.02 - 1.1 = 80mg
extracted w/ MeOH (Ag. knot w/ 60% DIPA in
meOH unnecessary) dried down from layer (yellow)
Bong - 60mg ($\frac{742 \text{ mg}}{\mu}$) = 35.5 mg peptide
Fmoc - 6 - 6 - F - L (50.14)

$$\text{Fmoc } 223.3 + 2(57.05) + 147.18 + 113.16 + 17 = 614.74 \frac{\text{mg}}{\text{mole}}$$

$$35.5 \text{ mg} / 614.74 \frac{\text{mg}}{\text{mole}} = 0.058 \text{ mmole} \times \left(\frac{2.35}{2.5} \right) = 5.5 \times 10^{-2} \text{ mmole DIPA}$$

$$(5.5 \times 10^{-2} \text{ mmole}) (379.25 \frac{\text{mg}}{\text{mole}}) = 20.9 \text{ mg HBTU}$$

$$(0.29 \text{ mmole}) (129.25 \frac{\text{mg}}{\text{mole}}) (1.4 \frac{\mu\text{L}}{\text{mg}}) = 52.5 \mu\text{L}$$

$$0.5 \text{ mL } (30 \text{ mM}) (614.74 \frac{\text{mg}}{\text{mole}}) = 9.22 \text{ mg } 9.5 \text{ mg peptide}$$

$$0.5 \text{ mL } (30 \text{ mM}) (379.25 \frac{\text{mg}}{\text{mole}}) = 5.68 \text{ mg HBTU}$$

$$0.5 \text{ mL } (60 \text{ mM}) (129.25 \frac{\text{mg}}{\text{mole}}) (0.742 \text{ mg}) = 5.6 \text{ mg DIPA}$$

Use 150 mM MeOH
9.15.04
Read and Understood By
13.16

Signed _____

Date _____

Signed _____

Signed _____

Date _____

16/29/04

Continued on Page 72

33.73 mg fmo₆-R₂H (2.54)

22.5 m, 173 TC (2.38 eV)

11.5 ~~pt~~ ~~1.024~~ (2.64 eq)

$$\frac{11.1}{2.69} = 4.12$$

71.7 ^{ser 7} 2 bird

4. Obwohl imf

Steel box @ 2:30 pm 09/5.04 50°C/150 Kpm

3:30 in " Kukul (3x pmf 2x moa) x

10 mm Run @ 150 Rpm $OD_{2.1} = 0.15$ To L₁ L₂

Winkel @ 50° ist zu groß für 10 mm - Prof 3x $\varnothing_{30} = 0,1$

inner RPR-run @ 150 RPM $OD_{300} = 2.9$

now need to get all the pipes off

Went @ 50°C 150 RPM Ruled 30 : left overnight

1.2 dm³ 2.1 dm³ 3x w/ perf $OD_{550} = 0.147$ Looks like

1st 10 min Break After Panel-Meeting

Complete peptide use $\frac{1}{3}$ of peptide solution. 27mg

Came out to 29mg & was clear so $\frac{1}{2}$ of the solution ✓

f_{11} 6 pin HSTU + 13 pc D.P.A. Generated 5 mV

Added to Sullivan @ 910 AM 7/16/07

2502 : 150 rpm

1:20 pm 9/16/14 Received & issued 4x w/ out

⑫

Revised @ 50% St Ag @ 125pm - 225pm

②

Rated 30 μ F 20mm Bank @ 50°C Rated 30 μ F

2

low R_{in} @ RT $OD_{301} = 0.11$

f) 10 mV 20^2 μ sec/mm @ RT 150 K/m $\phi_{201} = 1.75$

(R)

Printed 3x only

+ 19mg Naoc + 40mL Pipet 600 mL Perf

Continued on Page

50°C 150 rpm 3:30 9/14/17 Read and Understood By

Signed _____

Date

Signed

Date _____

3082.0


PROJECT

Removal from host @ 440 pm I placed in dish @ rt
(had been in the dish all along)

⑩ 920 pm Rinsed 3x in DMF ! left overnight
9/16/04

⑪ washed 3x DMF 10 min Rinse 0.35% =
600µl 60mM Semicarbazide + 20µl DAPI

Experimental Design 1st add Probe to patterned cells
This will cut any unwanted RNAs : block the antibody
on these patterned cells



1mJ add ^{quaternary} tertiary amine

730/2. 8mM FWHM $\lambda = 525\text{nm}$ 20x objective

1mJ for Y PC 2MJ for Both
0.1mJ for Y 1 Sec exposure

Probe Rxn $539.6 \frac{\text{mg}}{\text{mole}} = 337.36$ 21.1 mg Probe Rxn
33.73 mg

13 P	720
12 Y	660
10	660
10 Y	590
7 P	410
2 Y	420
7 Y	360
6 Y	310
5	240
4 Y	180
3 P	120
2 Y	60
1 P	0

⑫ 3x DMF, 10 min DMF, 3x DMF

+ 20µl Probe Rxn + 22.5µl H₂O + 21.7µl DAPI
50°C 150 rpm 12:07 pm 9.17.04 - 1:10 pm 9.17.04

⑬ Rinsed 3x DMF, 10 min @ 50°C 3x DMF, 10 min DMF

Continued on Page 73

Read and Understood By

[Signature] 9/27/04
Signed Date

[Signature] 10/29/04
Signed Date

OD_{571} of 10mM RuL = 0.086

① 10mM 20% reference $OD_{301} = 0.742$ (Scanned or Shunt)

② 300mM

$2.9 - 1.75 = 1.15$ not accessible to pftbe

have removed... my hair
bees on it and they want
recognition

$\frac{0.742}{2.9} = 25.82$ which is the number
of factors I detected 1 great

③ 600mM sucrose + 20mM RuL

1 sec @ 1mW to 0.1mW for 1st 6 rows (0.1mJ)

1 sec @ 1mW @ 1 sec (1mJ)

④ 300mM sucrose + 30mM RuL

Fractyl (or) 459.54 g/mole

$$\frac{539.6}{33.73} = \frac{459.54}{x}$$

+

28.73 mg Fractyl (or) + 22.5mg H3TU + 21.7mg RuL + 60mM sucrose

437 pm 50°C 150 RPM

⑤ 300mM

+ 600mM sucrose + 20mM RuL

1 sec @ 1mW (left side) w/ 2mJ (1 sec 2mW)
as before.

⑥ 300mM sucrose @ 50°C

Betaine monohydrate 94% (Sigma) $P_w = 135.2$

$$\frac{539.6}{33.7} = \frac{135.2}{x}$$

8.4mg Betaine
22.5mg H3TU
21.7mg RuL

Will get betaine to dissolve

300mM sucrose + 20mM RuL
at 50°C over the weekend

Read and Understood By

Continued on Page 75

Signed _____

Date

7/17/07

Signed _____

Signed

6/29/04

Date

PROJECT _____

Notebook No. _____

Continued From Page

74 75

9.20.04 Required 100% more to possible = 1.20 g
 Try 50% more, rather to be sure, it, went
 to minimum reaction error found.

① 500 μ l DMF + 80 μ l MeOH + 8.5 mg Acetone + 21.7 μ l D.P.A.

② + 22.5 μ l H₂O

3 ml added to chamber reacted @ RT 9.20.04 920 am 720

11:20 am 9.20.04 reacted 30 μ l DMF, 10 μ l DMF, 10 μ l DMF

OP₃₀₁ = 0.086

④ 10 ml 20% p.p.s. solution

OP₃₀₁ = 0.39

⑫ mixed 30 DMF 15 ml DMF 3X DMF

~ 50% of
 ~ 2 ml

Antibody hybridized

3X wash w/ PBS Buffer pH = 7.2 30 ml

100 μ l of L-8516 UV-crosslinked Antibody (1235.4)

+ 900 μ l PBS

2:11 pm 9.20.04 - 3:52 pm

TM 9.20.04

Reacted 3X with 0.05% Tween 20 / PBS

4 pm added 1/1000 Dilution of Antisense

5 pm mixed 20 μ l 0.05% Tween 20 in PBS

1 ml of 488 nm excitation 75% SL

14 ml got the opposite effect I was expect.

Could be that the antibody binds YP66FL or it
 more likely is from the Fluorescent photo product.

Conclusion | The most interesting part so far is that E can
 getting ~ 100% yield @ (2mJ) I saw much more

Fluorescence which goes in getting close to

100% w/ (2mJ)

Next Read and Understood for more ms

Continued on Page 76

Signed

Date

Signed

Date

Scraped off the two tubes of the vial pattern
 & placed the polymerized treatments into separate microfuge
 tubes added 200 μ l 50:50 TFA in Dean w/ 3% Triisopropylsilane

10:55 am 9.21.04

1:30 pm still yellow still around without? don't know


Observation the fibers are yellow, both for the treatment &
 control. I hope they will become clear after the cleanse?

did not change in 1:1 Dean: TFA

1:35 pm added neat TFA & immediately the
 fibers went from yellow after small
 amounts added 0.3 μ l Triisopropylsilane

not TFA, polymer floats in TFA so it
 was dispersed so it GLIA see the color
 added 100 μ l acetonitrile & color less color

but brown & green to one for m.s

 9/22/04

1 tan 9/23/04

2 Repeated cleanse on 2 rows of B-66 PL mass
 in one row NVOG-66 PL - mass

in 12 TFA supernatant silane in TFA


cleaned lower array (466 PL & 4766 PL) in 1:1

TFA = Acetonitrile w/ 0.5% Triisopropylsilane same
 when polymer may have fallen apart

MALDI-MS shows m/z = 440.32 peak which is 1.00 off from predicted fragment
 structure control shows m/z = 653.259. NVOG 66 PL - mass
 done 10/29/04

Continued on Page

Read and Understood By



Signed

9.22.04


Date





Signed

10/29/04


Date




Calibrated MALDI-TOF MS spectrum showing ions formed from photopatterned N-Tris(2,4,6-trimethoxyphenyl)phosphine-GGFL ($m/z=964.4$ Da) peptide. Inset image of photopatterned array of Texas Red sulfonylethylphosphine (Red) and fluorescein isothiocyanate (green). N-Tris(2,4,6-trimethoxyphenyl)phosphine (TMPP) is facilitates product detection and formation of a ions for post source decay analysis. Analytical Biochemistry 268, 305-317 (1999).





Calibrated MALDI-TOF MS spectrum of observed isotopic distribution for the $m/z=964.4$ Da ion vs. those predicted for the TMPP-GGFL [C48H63N5O14P] (bottom).




Uncalibrated MALDI-TOF MS post source decay showing the a1 (TMPP-G), a2 (TMPP-GG), a3 (TMPP-GGF), and primary ion m/z=964.8 Da of the TMPP-GGFL peptide.




Calibrated MALDI-TOF MS spectrum showing ions formed from control (not irradiated) areas.



sample B, betaine-GGFL-amido
 D:\1027\ref_0004.dat
 Acquired: 11:15:00, September 28, 2004





**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☒ **BLACK BORDERS**

☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**

☐ **FADED TEXT OR DRAWING**

☒ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**

☐ **SKEWED/SLANTED IMAGES**

☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**

☐ **GRAY SCALE DOCUMENTS**

☒ **LINES OR MARKS ON ORIGINAL DOCUMENT**

☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**

☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.

From the INTERNATIONAL BUREAU

PCTNOTIFICATION CONCERNING
SUBMISSION OR TRANSMITTAL
OF PRIORITY DOCUMENT

(PCT Administrative Instructions, Section 411)

Date of mailing (day/month/year) 19 July 2005 (19.07.2005)		IMPORTANT NOTIFICATION
Applicant's or agent's file reference 112624.00138 PCT		
International application No. PCT/US2005/015764	International filing date (day/month/year) 06 May 2005 (06.05.2005)	
International publication date (day/month/year)	Priority date (day/month/year) 06 May 2004 (06.05.2004)	
Applicant ARIZONA BOARD OF REGENTS, acting for and on behalf of, Arizona State University et al		

1. By means of this Form, which replaces any previously issued notification concerning submission or transmittal of priority documents, the applicant is hereby notified of the date of receipt by the International Bureau of the priority document(s) relating to all earlier application(s) whose priority is claimed. Unless otherwise indicated by the letters "NR", in the right-hand column or by an asterisk appearing next to a date of receipt, the priority document concerned was submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b).
2. *(If applicable)* The letters "NR" appearing in the right-hand column denote a priority document which, on the date of mailing of this Form, had not yet been received by the International Bureau under Rule 17.1(a) or (b). Where, under Rule 17.1(a), the priority document must be submitted by the applicant to the receiving Office or the International Bureau, but the applicant fails to submit the priority document within the applicable time limit under that Rule, the attention of the applicant is directed to Rule 17.1(c) which provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.
3. *(If applicable)* An asterisk (*) appearing next to a date of receipt, in the right-hand column, denotes a priority document submitted or transmitted to the International Bureau but not in compliance with Rule 17.1(a) or (b) (the priority document was received after the time limit prescribed in Rule 17.1(a) or the request to prepare and transmit the priority document was submitted to the receiving Office after the applicable time limit under Rule 17.1(b)). Even though the priority document was not furnished in compliance with Rule 17.1(a) or (b), the International Bureau will nevertheless transmit a copy of the document to the designated Offices, for their consideration. In case such a copy is not accepted by the designated Office as the priority document, Rule 17.1(c) provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.

<u>Priority date</u>	<u>Priority application No.</u>	<u>Country or regional Office or PCT receiving Office</u>	<u>Date of receipt of priority document</u>
06 May 2004 (06.05.2004)	60/569,370	US	20 June 2005 (20.06.2005)
10 September 2004 (10.09.2004)	60/608,774	US	20 June 2005 (20.06.2005)
29 October 2004 (29.10.2004)	60/623,181	US	20 June 2005 (20.06.2005)

The International Bureau of WIPO
34, chemin des Colombettes
1211 Geneva 20, Switzerland

Facsimile No. +41 22 338 82 70

Form PCT/IB/304 (January 2004)

Authorized officer

Paulette BOCCARD

Facsimile No. (41-22) 338.87.40

Telephone No. +41 22 338 8147

CIP17012

From the INTERNATIONAL BUREAU

PCTNOTIFICATION CONCERNING
SUBMISSION OR TRANSMITTAL
OF PRIORITY DOCUMENT

(PCT Administrative Instructions, Section 411)

Date of mailing (day/month/year) 19 July 2005 (19.07.2005)		To: ATKINS, Robert, D. Quarles & Brady Streich Lang, LLP One Renaissance Square Two North Central Avenue Phoenix, AZ 85004 ETATS-UNIS D'AMERIQUE	
Applicant's or agent's file reference 112624.00138 PCT			IMPORTANT NOTIFICATION
International application No. PCT/US2005/015764	International filing date (day/month/year) 06 May 2005 (06.05.2005)		
International publication date (day/month/year)	Priority date (day/month/year) 06 May 2004 (06.05.2004)		
Applicant ARIZONA BOARD OF REGENTS, acting for and on behalf of, Arizona State University et al			

1. By means of this Form, which replaces any previously issued notification concerning submission or transmittal of priority documents, the applicant is hereby notified of the date of receipt by the International Bureau of the priority document(s) relating to all earlier application(s) whose priority is claimed. Unless otherwise indicated by the letters "NR", in the right-hand column or by an asterisk appearing next to a date of receipt, the priority document concerned was submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b).
2. (If applicable) The letters "NR" appearing in the right-hand column denote a priority document which, on the date of mailing of this Form, had not yet been received by the International Bureau under Rule 17.1(a) or (b). Where, under Rule 17.1(a), the priority document must be submitted by the applicant to the receiving Office or the International Bureau, but the applicant fails to submit the priority document within the applicable time limit under that Rule, the attention of the applicant is directed to Rule 17.1(c) which provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.
3. (If applicable) An asterisk (*) appearing next to a date of receipt, in the right-hand column, denotes a priority document submitted or transmitted to the International Bureau but not in compliance with Rule 17.1(a) or (b) (the priority document was received after the time limit prescribed in Rule 17.1(a) or the request to prepare and transmit the priority document was submitted to the receiving Office after the applicable time limit under Rule 17.1(b)). Even though the priority document was not furnished in compliance with Rule 17.1(a) or (b), the International Bureau will nevertheless transmit a copy of the document to the designated Offices, for their consideration. In case such a copy is not accepted by the designated Office as the priority document, Rule 17.1(c) provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.

<u>Priority date</u>	<u>Priority application No.</u>	<u>Country or regional Office or PCT receiving Office</u>	<u>Date of receipt of priority document</u>
06 May 2004 (06.05.2004)	60/569,370	US	20 June 2005 (20.06.2005)
10 September 2004 (10.09.2004)	60/608,774	US	20 June 2005 (20.06.2005)

The International Bureau of WIPO 34, chemin des Colombettes 1211 Geneva 20, Switzerland Facsimile No. +41 22 338 82 70	Authorized officer Paulette BOCCARD Facsimile No. (41-22) 338.87.40 Telephone No. +41 22 338 8147
---	---

From the INTERNATIONAL BUREAU

PCTNOTIFICATION CONCERNING
SUBMISSION OR TRANSMITTAL
OF PRIORITY DOCUMENT

(PCT Administrative Instructions, Section 411)

Date of mailing (day/month/year) 19 July 2005 (19.07.2005)	IMPORTANT NOTIFICATION
Applicant's or agent's file reference 112624.00138 PCT	
International application No. PCT/US2005/015764	
International publication date (day/month/year)	
Applicant ARIZONA BOARD OF REGENTS, acting for and on behalf of, Arizona State University et al	

To:

ATKINS, Robert, D.
Quarles & Brady Streich Lang, LLP
One Renaissance Square
Two North Central Avenue
Phoenix, AZ 85004
ETATS-UNIS D'AMERIQUE

- By means of this Form, which replaces any previously issued notification concerning submission or transmittal of priority documents, the applicant is hereby notified of the date of receipt by the International Bureau of the priority document(s) relating to all earlier application(s) whose priority is claimed. Unless otherwise indicated by the letters "NR", in the right-hand column or by an asterisk appearing next to a date of receipt, the priority document concerned was submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b).
- (If applicable)* The letters "NR" appearing in the right-hand column denote a priority document which, on the date of mailing of this Form, had not yet been received by the International Bureau under Rule 17.1(a) or (b). Where, under Rule 17.1(a), the priority document must be submitted by the applicant to the receiving Office or the International Bureau, but the applicant fails to submit the priority document within the applicable time limit under that Rule, the attention of the applicant is directed to Rule 17.1(c) which provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.
- (If applicable)* An asterisk (*) appearing next to a date of receipt, in the right-hand column, denotes a priority document submitted or transmitted to the International Bureau but not in compliance with Rule 17.1(a) or (b) (the priority document was received after the time limit prescribed in Rule 17.1(a) or the request to prepare and transmit the priority document was submitted to the receiving Office after the applicable time limit under Rule 17.1(b)). Even though the priority document was not furnished in compliance with Rule 17.1(a) or (b), the International Bureau will nevertheless transmit a copy of the document to the designated Offices, for their consideration. In case such a copy is not accepted by the designated Office as the priority document, Rule 17.1(c) provides that no designated Office may disregard the priority claim concerned before giving the applicant an opportunity, upon entry into the national phase, to furnish the priority document within a time limit which is reasonable under the circumstances.

<u>Priority date</u>	<u>Priority application No.</u>	<u>Country or regional Office or PCT receiving Office</u>	<u>Date of receipt of priority document</u>
06 May 2004 (06.05.2004)	60/569,370	US	20 June 2005 (20.06.2005)

The International Bureau of WIPO
34, chemin des Colombettes
1211 Geneva 20, Switzerland

Facsimile No. +41 22 338 82 70

Authorized officer

Paulette BOCCARD

Facsimile No. (41-22) 338.87.40
Telephone No. +41 22 338 8147